

Supplementary data on method for analysis of
2-Aminoadipic acid (2AAA).

Method based on article

Midttun et al (2016), PMID 27715010.

Material

2-Aminoadipic acid (purity 98%) was obtained from Fluorchem Ltd,
Unit 15, Graphite Way Hadfield Glossop Derbys, SK13 1QH UK.
2-Aminoadipic acid-d3 (purity >98%) was obtained from C/D/N
Isotopes Inc., 88 Leacock Street Pointe-Claire, Quebec, Canada.

Instrumentation

Agilent 7010B GC/TQ and Agilent 8890 GC System.

Chromatography and detection

GC-MS/MS; positive-ion multiple reaction monitoring (MRM);
retention time = 4.31 min.

2-Aminoadipic acid precursor ion = 202.0 m/z; product ion =
114.0 m/z.

2-Aminoadipic acid-d3 precursor ion = 205.0 m/z; product ion =
115.0 m/z.

Method performance

Linear range: 0.2 - 100 $\mu\text{mol/L}$.

Linearity: r^2 : 0.997.

LOD (S/N >5): 0.2 $\mu\text{mol/L}$.

Within-day CV: 3-5 %.

Between-day CV: 5-7 %.