

Supplementary data on method for analysis of
3-Hydroxybutyrate (bHB).

Method based on article

Midttun et al (2016), PMID 27715010.

Material

3-Hydroxybutyrate (purity 95%) was obtained from
Sigma-Aldrich, St.Louis, MO 63103 USA or 89555 Steinheim
Germany.

3-Hydroxybutyrate-d2 (purity 99%) was obtained from C/D/N Isotopes
Inc., 88 Leacock Street Pointe-Claire, Quebec, Canada.

Instrumentation

Agilent 7010B GC/TQ and Agilent 8890 GC System.

Chromatography and detection

GC-MS/MS; positive-ion multiple reaction monitoring (MRM);
retention time = 2.95 min.

3-Hydroxybutyrate precursor ion = 145.0 m/z; product ion = 69.0 m/z.

3-Hydroxybutyrate-d2 precursor ion = 147.0 m/z; product ion =
70.0 m/z.

Method performance

Linear range: 1 - 1000 $\mu\text{mol/L}$.

Linearity: r^2 : 0.998.

LOD (S/N >5): 1 $\mu\text{mol/L}$.

Within-day CV: 2-4 %.

Between-day CV: 3-5 %.