

Supplementary data on method for analysis of 3-Hydroxyisobutyrate (3-HIB).

Method based on article

Midttun et al (2016), PMID 27715010.

Material

3-Hydroxyisobutyrate (purity $\geq 96\%$) was obtained from Sigma-Aldrich, St.Louis, MO 63103 USA or 89555 Steinheim Germany.

3-Hydroxyisobutyrate-d5 (purity 99+%) was obtained from BUCHEM BW, Minden 60/62 7327 AW Apeldoorn, The Netherlands.

Instrumentation

Agilent 7010B GC/TQ and Agilent 8890 GC System.

Chromatography and detection

GC-MS/MS; positive-ion multiple reaction monitoring (MRM);
retention time = 2.98 min.

3-Hydroxyisobutyrate precursor ion = 145.0 m/z; product ion =
69.0 m/z.

3-Hydroxyisobutyrate-d5 precursor ion = 150.0 m/z; product ion =
74.0 m/z.

Method performance

Linear range: 1 - 1000 $\mu\text{mol/L}$.

Linearity: r^2 : 0.998.

LOD (S/N >5): 1 $\mu\text{mol/L}$.

Within-day CV: 2-4 %.

Between-day CV: 3-5 %.