

Supplementary data on method for analysis of Carboxymethyllysine (CML).

Method based on article

Middtun et al (2016), PMID 27715010.

Material

Carboxymethyllysine (purity $\geq 98\%$) was obtained from Carbosynth Limited, 8&9 Old Station Business Park, Compton, Berkshire, RG20 6NE, UK.

Carboxymethyllysine-d3 (purity $>98\%$) was obtained from C/D/N Isotopes Inc., 88 Leacock Street Pointe-Claire, Quebec, Canada.

Instrumentation

Agilent 7010B GC/TQ and Agilent 8890 GC System.

Chromatography and detection

GC-MS/MS; positive-ion multiple reaction monitoring (MRM); retention time = 6.93 min.

Carboxymethyllysine precursor ion = 227.0 m/z; product ion = 198.0 m/z.

Carboxymethyllysine-d3 precursor ion = 229.0 m/z; product ion = 200.0 m/z.

Method performance

Linear range: 0.03 - 100 $\mu\text{mol/L}$.

Linearity: r^2 : 0.999.

LOD (S/N >5): 0.03 $\mu\text{mol/L}$.

Within-day CV: 3-5 %.

Between-day CV: 3-5 %.