

## ALFA elution peptide

The ALFA elution peptide can be used for competitive elution of ALFA-tagged proteins from magnetic or non-magnetic ALFA Selector<sup>PE</sup> and ALFA Selector<sup>CE</sup> resins.

The peptide sequence is Acetyl-PSRLEEELRRRLTEP-Amide.

The protocol below refers to the following product:

Cat. No.	Product Name	Lyophilized from	Reconstitution	Final conc.
N1520	ALFA elution peptide	ddH <sub>2</sub> O	250 µL	40 mg/mL

The ALFA elution peptide is shipped as lyophilized powder at ambient temperature. The lyophilized reagent can be stored at 2-8°C for up to 12 months.

## Protocol: Reconstitution and Storage of ALFA elution peptide

- 1. Prepare sterile deionized water. If applicable, we recommend including 0.1% sodium azide as a preservative.
- 2. Open the vial containing the lyophilized ALFA elution peptide.
- 3. Add 250 µL of sterile deionized water.
- Vortex or mix at room temperature for 5-10 min.
  If the material is not completely dissolved, it may help to increase the temperature to 37°C an mix for up to 30 min.
- 5. Briefly spin down the vial for 2 min at 1000 xg.
- 6. Distribute into aliquots. Use small tubes and avoid aliquots below 20  $\mu L.$
- 7. Storage after reconstitution:
  - Aliquots may be stored at -20°C or below.
- Before usage for competitive elution from ALFA Selector resins, the concentrated peptide solution obtained after reconstitution (40 mg/mL corresponding to 20 mM peptide) has to be diluted in a buffered solution. Depending on the specific application, the final peptide concentration may vary from 200-1000 μM.

## Only for research applications, not for diagnostic or therapeutic use!