

Immunostaining of cultured cells using FluoTag[®] reagents

Materials (not included):

- Phosphate Buffered Saline (PBS) pH 7.4
- Paraformaldehyde (PFA):
4% PFA in PBS pH 7.4, freshly prepared ⁽¹⁾
- Quenching Solution (QS):
0.1 M Glycine in PBS pH 7.4
or: 0.1 M NH₄Cl in PBS pH 7.4
- Blocking & Permeabilization buffer (BPB) ⁽²⁾:
10% Normal Goat Serum (NGS) + 0.1% Triton X-100 in PBS
or: 2% Bovine Serum Albumin (BSA) + 0.1% Triton X-100 in PBS
- FluoTag[®] Dilution Buffer (FDB) ⁽²⁾:
3% NGS + 0.1% Triton X-100 in PBS
or: 1% BSA + 0.05% Triton X-100 in PBS
- High-salt PBS:
PBS supplemented with 0.5 M NaCl

Note: Below are two versions of the same IF protocol. The "Bench Protocol" is a concise version, while the "Detailed Protocol" provides an extended explanation.

Bench Protocol: IF of cultured cells using FluoTag[®] reagents

1. Fix cells with 4% PFA for 30 min at RT.
FluoTag[®] reagents are also compatible with methanol fixation. In this case, step 2 can be omitted.
2. Quench with PBS supplemented with 0.1 M glycine or 0.1 M NH₄Cl for 10 min at RT.
3. Wash once with PBS.
4. Permeabilize and block for 15 min with 10% Natural Goat Serum (NGS) and 0.1% Triton X-100 in PBS.
5. Dilute FluoTag[®] with 3% NGS and 0.1% Triton X-100 in PBS.
To obtain optimal results for different target proteins and expression levels, the dilution of FluoTag[®] products might need to be adjusted. The recommended dilution is thus only a starting point for optimization.
6. Incubate fixed cells with the diluted FluoTag[®] for 60 min at RT.
7. Wash 3 times for 5 min each with 1 mL of PBS.
8. Optional: Wash once with high-salt PBS (PBS + 500 mM NaCl) followed by PBS.
9. Shortly dip coverslip in water before mounting.
We recommend using Mowiol as a mounting medium.

Detailed Protocol: IF of cultured cells using FluoTag[®] reagents

- Note:
- The example below is based on a 12-well plate.
 - Please adapt the protocol to your experimental conditions.

1. Wash cells gently using PBS (e.g. 1 mL of PBS per well).
2. Add 1 mL of 4% PFA per well and incubate at room temperature (30 min, RT).
3. Remove PFA and dispose according your laboratory rules.
4. Briefly rinse with 1 mL QS per well.
5. Add 1 mL of fresh QS per well and shake gently on an orbital shaker (10 min, RT).
6. Remove QS.
7. Briefly rinse with 1 mL of PBS per well.
8. Add 1 mL of BPB per well and shake gently (15 min, RT).
9. During this time prepare the FluoTag[®] working solution. Make sure to prepare sufficient volume for all reactions (e.g. 5 mL for a full 12 well plate).
 - a. Vortex FluoTag[®] stock solution shortly and centrifuge for 2 min at 10.000 *xg*.
 - b. Dilute the FluoTag[®] reagent in FluoTag[®] Dilution Buffer ⁽³⁾.
10. Remove BPB solution from wells.
11. Add 400 μ L per well of the FluoTag[®] working solution. Incubate for 60 min. with gentle shaking at RT and protected from light.
12. Remove the FluoTag[®] working solution from well.
13. Rinse once with 1 mL of PBS per well.
14. Wash with 1 mL of PBS per well and shake the plate gently for 5 min at RT and protected from light.
15. Repeat the previous step 2 times.
16. Optional step:
 - a. Wash once for 5 min with high-salt PBS.
 - b. Briefly rinse with PBS.
17. Before mounting, rinse once with water to remove the excess of salt.

Remarks

- (1) FluoTag[®] products are also compatible with methanol fixation. Fixation protocols using glutaraldehyde are not recommended.
- (2) We recommend using blocking and FluoTag[®] Dilution Buffers prepared with Normal Goat Serum (NGS).
- (3) To obtain optimal results for different target proteins and expression levels, the dilution factor might need to be adjusted. The recommended dilution specified in the data sheet is thus only a starting point for further optimizations.

For further information please contact us at info@nano-tag.com