



The Lanluma non-surgical *buttock lift*

The revolutionary collagen
stimulating injectable.

||||| LANLUMA

Get results comparable
to buttock surgery with
minimal risk and downtime.

What's inside

05

Introducing
Lanluma

09

The procedure
at a glance

06

Curvy buttocks
you'll want to
show off

10

How does
Lanluma work?

08

What happens
during a Lanluma
treatment?



Introducing Lanluma

Lanluma is the collagen stimulating injectable that smoothes, reshapes and sculpts beautiful buttocks.

This innovative filler reshapes and activates natural collagen production, creating volume that enhances your body's curves and contours.

PLLA particles get to work deep within the skin to activate your body's own natural collagen response – so the volume it creates is long-lasting. Injectable PLLA, Lanluma's main ingredient, provides results that appear gradually and last more than two years.* Lanluma is the perfect choice to boost your appearance, and as a result, your confidence.

Did you know?

Lanluma's main ingredient is poly-L-lactic acid (PLLA). This clever component has been used in medicine for over 30 years for everything from dissolvable stitches to treating scarring.

The lowdown on the Lanluma buttock boost

- Lanluma provides comparable results to buttock surgery with minimal risk and downtime.
- This amazing injectable activates your own natural collagen deep within your skin to create volume where you want it most.
- Suitable for adults of all ages and gender identities.
- You get beautifully natural-looking results. Injectable PLLA, Lanluma's main ingredient, provides results that last more than two years.*
- It promises a gradual, continuous improvement in both volume and skin texture as your treatment plan progresses and your natural collagen production is stimulated.

*J Am Acad Dermatol 2010;62:448-62

Curvy buttocks you'll want to show off

Lanluma enhances the body's curves and contours, and has been specially designed to help you get the full, smooth, peachy buttocks you've dreamed of.

What's more, the treatment is suitable for adults of all gender identities and ages.

How does the Lanluma buttock boost work?

This innovative, minimally invasive treatment sets off a deep tissue regeneration reaction. It stimulates collagen, adding volume, improving skin texture and recontouring your buttocks for a more shapely silhouette.

Lanluma helps to correct and enhance in many ways:


Enhances curves in flat or saggy buttocks.

Improves skin tightness and texture on the buttocks.

Can help to reduce cellulite and dimples.

Improves the shape and volume of the buttocks, whatever your age or gender.

Corrects natural and post-surgical concerns such as atrophy and lipoatrophy.



“ Lanluma stimulates collagen, adding volume, improving skin texture and recontouring your buttocks for a more shapely silhouette. ”



What happens during a Lanluma treatment?

It all starts with a conversation. Before any treatment begins, your clinician will encourage you to talk through what you'd like to achieve.

They'll look at the areas you'd like to transform, and discuss what's possible.

Then, if you are happy to go ahead, treatment can proceed. Unlike other collagen stimulating injectables – for which you will have to schedule another session – Lanluma is quick to prepare so it can be injected from the first consultation.

After making you comfortable, they will sterilise the area to be treated, and then apply a local anaesthetic. This makes the treatment virtually painless. They will carefully inject Lanluma where you would like to create volume – and that's it.

You'll feel fine right after the procedure, though there may be a little redness and swelling for a few hours.

Depending on the area treated and how your body reacts, creating the volume you're after will take 2-3 sessions, and results will show in the months afterwards as your natural collagen develops.

Aftercare

Downtime is minimal, and your physician will explain fully and provide an easy aftercare guide to follow. You'll need to keep up a careful aftercare routine of regular massage for best results – 10 minutes, 2 x a day for two weeks after your last session. It's easy to maintain though – most people set a regular phone alarm, and some even find it a relaxing way to take some 'me time' in the day. You will also need to avoid swimming, sauna and vigorous exercise for a week after treatment, and sunbathing (UV) for 10 days.

The procedure at a glance



It's an outpatient treatment that can be typically performed in **under an hour**.



It's virtually painless (performed under local anaesthetic) with minimal downtime – patients can resume normal activity **same-day**.



It takes between **2-3** treatment sessions, with results appearing gradually over **2-3** months.



Aftercare is easy – massage for 10 minutes, **2 x a day for two weeks**.



How does Lanluma work?

Lanluma is an injectable poly-L-lactic acid (PLLA) collagen stimulant.

It can be used to activate natural collagen in two ways: either to increase fullness and add volume; or to create a more youthful appearance by helping the body naturally correct sagging skin, scarring and dimples.

Why is it different?

Lanluma isn't the same as HA fillers, because it doesn't simply 'fill up' the area where you would like more volume. In fact, PLLA particles get to work deep within the skin to activate the body's own natural collagen response – so the volume created is longer-lasting. It also improves your skin texture from within.

Safer than surgery

Lanluma differs from buttock lift surgery because it's very safe and minimally invasive. PLLA has been used in medicine for over 30 years. There's no risk of general anaesthetic, and chances of complications such as infection are very low compared to buttock surgery.

Minimal downtime

Finally, a Lanluma buttock boost demands virtually no downtime. Though you might have a little redness and swelling after each session, you won't need to take time off work. You'll need to avoid sunbathing, gym and sauna for a few days after each session and keep up your massage routine, but that's all – simple, safe and effective, you'll soon be wearing your new curves and confidence with pride.

*“ Simple, safe and effective,
you’ll soon be wearing your new
curves and confidence with pride.”*



||||| LANLUMA

SINCLAIR PHARMA LTD.

1ST FLOOR
WHITFIELD COURT
30-32 WHITFIELD STREET
LONDON
W1T 2RQ

www.lanlumafiller.com

SINCLAIR

The Sinclair logo features the word "SINCLAIR" in a serif font, positioned to the left of a stylized, dark blue graphic element. This graphic element consists of two curved lines that form a partial circle or a stylized letter 'S', with the top curve extending upwards and the bottom curve extending downwards.