BORT Generation

Kniebandage / Kniebandage Plus

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Generation

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
REF 114 490 G/215 420

Illustrationen
Figures

1

2

3

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	04
EN	english	Instructions for use	08
FR	français	Mode d'emploi	12
ES	español	Instrucciones de uso	16
IT	italiano	Istruzioni per l'uso	20
NL	nederlands	Gebruiksaanwijzing	24
CS	český	Návod k použití	28
ET	eesti	Kasutusjuhend	32
PL	polski	Instrukcja użytkowania	36
RO	românesc	Instrucțiuni de utilizare	40

 PDF: ga.bort.com

BORT Generation Kniebandage **BORT Generation Kniebandage Plus**

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BORT Generation Knee Support **BORT Generation Knee Support Plus**

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a knee support for compression of soft tissue with a local pressure pad around the kneecap consisting of flexible material as well as a circumferential sealing tape for improved stabilisation on the thigh as an additional functional element. (additional sealing tape only for Generation Knee Support Plus REF 114 490 G).

Indications

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Before fitting the latch plates on the thigh and for REF 114 490 G, open the additional sealing tape.

1 Open the velour fastening on the knee support thigh. The Generation Knee Support Plus (REF 114 490 G) has an additional circumferential harness, please also open it. If required, the harness can be unthreaded. The joint straps should always be located on the outer side of the leg.

2 Pull the knee support as far as possible over the knee so that the silicone pad encloses the kneecap.

3 The Generation Knee Support Plus (REF 114 490 G) additionally has a circumferential harness, please guide it through the guide loops and fix it. The joint straps should always be located on the outer side of the leg. Finally, check the fit and correct positioning of the knee support, so that the support supports and compresses but doesn't constrict.

Removal

To remove the support, open the sealing tape (only for the Generation Knee Support Plus REF 114 490 G) and the latch plates on the thigh, grasp the lower edge of the support at the side in the vicinity of the splints and pull it down over the foot.

Material composition

Polyamide (PA), elastodiene/latex (LA), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash Do not clean chemically Do not bleach

Do not dry in a tumble dryer Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 01.2020

Medical device |  Single patient – multiple use

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