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
REF 121 980

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D121980I2020-05I004 ML

DE BORT OmoSAT

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EN BORT OmoSAT

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined position.

Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of distortion, contusion, subcapital humerus fracture, signs of paralysis

Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/important notes ⚠

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on our around it as they can destroy the material

Fitting instructions

- Place the forearm in the arm pouch and position your elbow in the semi-circular stitched down area.
- Position the training ball in the metacarpus area and fix it using the Velcro area in the arm pouch.
- The arm pouch is closed at the height of the central lower arm region using the short forearm strap included. This strap fixes the forearm in the pouch.
- The shoulder strap runs from the elbow area ring over the opposite side of the shoulder height through the lower arm area ring. The strap system Velcro fasteners enable individual setting of the desired resting position.

5 The torso strap is used for additional fixing to the body. It runs from the elbow area ring horizontally around the back through the lower arm area ring. Its circumference can also be set individually using the Velcro fasteners.

Removal

To remove the OmoSAT shoulder joint brace, open all Velcro fasteners, loosen all Velcro fasteners in the reverse order, then remove the brace.

Material composition

Polyamide (PA), PUR foam (PUR), polyester (PES), cotton (CO)
The sewn in textile label provides the precise material composition.

Cleaning information

☑ Delicate wash ☒ Do not clean chemically ☒ Do not bleach
☒ Do not dry in a tumble dryer ☒ Do not iron
Close the Velcro fastening to avoid damage to other laundry items.
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2019

Medical device | Single patient – multiple use

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