BORT Generation Knieorthese

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Generation

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

REF 215 400

Illustrationen
Figures


1


2

3

4

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	04
EN	english	Instructions for use	08
FR	français	Mode d'emploi	12
ES	español	Instrucciones de uso	16
IT	italiano	Istruzioni per l'uso	21
NL	nederlands	Gebruiksaanwijzing	26
CS	český	Návod k použití	30
ET	eesti	Kasutusjuhend	34
PL	polski	Instrukcja użytkowania	38
RO	românesc	Instrucțiuni de utilizare	42

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a support for supporting the physiological guidance of the knee joint with circumferential Velcro fastener restraints and side guidance splints with joints which are adjustable in their scope of movement.

Indications

Degenerative diseases, e.g. mild to moderate lateral ligament instability of the knee joint, gonarthrosis, meniscus injuries.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Before putting the support on, the side joint splints must be adapted to the contour of the leg by an expert. For forming, remove the splints from the textile and align to the lateral course of the patient's leg by bending it up (if applicable using a bending iron). Then, set the extension and flexion restrictions on the medial and lateral toothed segment articulation. The extension is factory set to 0° , flexion is unrestricted. A different extension stop can be attached by loosening the assembly screw, removing the limit stopper inserted ex-works, inserting the desired limit stopper and attaching it with the assembly screw. You can work accordingly with flexion wedges.

- 1** For fitting, leave the two Velcro straps and the BORT Generation knee brace padding supports below them open.
- 2** When sitting, lay the BORT Generation knee brace on the stretched knee joint so that the patella recess is directly placed on the kneecap.
- 3** Then guide the padding supports around the thigh and shank and close them.
- 4** Finally guide the two Velcro straps through the guide loops and close them with the desired tension.

Removal

To remove the aid, open the circumferential Velcro straps and then the surface padding supports. Then, remove the support from the leg. Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

Material composition

Polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

Cleaning information

Before washing, remove the joint splints and close the Velcro fasteners in order to avoid damaging other items of clothing.

 Delicate wash  Do not clean chemically  Do not bleach

 Do not dry in a tumble dryer  Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 12.2019

Medical device |  Single patient – multiple use

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