# BORT Generation Rückenbandage

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Generation

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**REF** 215 850

## Sprachen Languages

<b>DE</b>	deutsch	Gebrauchsanweisung .....	03
<b>EN</b>	english	Instructions for use .....	07
<b>FR</b>	français	Mode d'emploi .....	11
<b>ES</b>	español	Instrucciones de uso .....	15
<b>IT</b>	italiano	Istruzioni per l'uso .....	19
<b>NL</b>	nederlands	Gebruiksaanwijzing .....	23
<b>CS</b>	český	Návod k použití .....	27
<b>ET</b>	eesti	Kasutusjuhend .....	31
<b>PL</b>	polski	Instrukcja użytkowania .....	35
<b>RO</b>	românesc	Instrucțiuni de utilizare .....	39

 PDF: [ga.bort.com](http://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

## **Intended purpose**

This medical device is a brace with a back pad and tension elements for stabilising the lumbar spine.

## **Indications**

Degenerative lumbar changes, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenosis, radicular irritation syndrome, condition following intervertebral disc surgery.

## **Contraindications**

Impaired lymphatic drainage including unclear soft tissue swelling, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

## **Application risks/Important notes**

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds

- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

### **Fitting instructions**

When you receive your BORT Generation Back Support, the front fastener is closed. Open it fully. The pad which is split into two in the back area must be positioned with the respective wider part of the pad in the direction of the head. Place both pads respectively on the outer edge of the support cover. Place the back support around the abdomen and close the BORT Generation Back Support applying very little force, using the attached hand straps. In order to close the BORT Generation Back Support fully, slip into the hand straps and pull them forwards simultaneously and evenly, applying very little force. Following this, attach the hand straps to the front fastener or the velour surfaces at the side using the micro Velcro provided. Observe here that the support is not closed too tightly. Within the course of this fitting process, the BORT Generation Back Support is fully fitted to the anatomy of the body and interlocks. If the lower part of the support ends with the hip-bend when sitting, the support is correctly positioned.

### **Removal**

To remove the support, loosen the pulley hand straps and then open the front fastener.

### **Material composition**

Polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

## Cleaning information

☞ Hand-washing ☒ Do not clean chemically ☒ Do not bleach

☒ Do not dry in a tumble dryer ☒ Do not iron

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:  
[www.bort.com/conformity](http://www.bort.com/conformity)

Status: 11.2019

Medical device |  Single patient – multiple use

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