BORT EpiBasic Sport

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


Sport

Illustrationen
Figures

1

2

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	04
EN	english	Instructions for use	08
FR	français	Mode d'emploi	11
ES	español	Instrucciones de uso	15
IT	italiano	Istruzioni per l'uso	19
NL	nederlands	Gebruiksaanwijzing	23
CS	český	Návod k použití	27
ET	eesti	Kasutusjuhend	30
PL	polski	Instrukcja użytkowania	34
RO	românesc	Instrucțiunile de utilizare	38

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a circumferential compression elbow compression support with elastic material fabric pads and elastic profile pads.

Indications

Chronic post-traumatic or post-operative soft tissue irritation in the elbow joint area, epicondylitis, tendomyopathy.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Place the support on the elbow affected and position the knitted comfort zone in the antecubital fossa at the same time. The support is correctly positioned if you can feel the outer and inner condyle through the hollows in the silicone insert.

1 Stabilo strap with the easyClick system for individual dosage of the pad pressure: **1** and **3** attachment module, **2** supporting module, **4** adjustment module. Adapt the pad pressure dosage individually by moving the adjustment module (**4**).

2 When you receive your support, the attachment Stabilo strap module (**1**) in the supporting module (**2**) is attached on the support. Guide the Stabilo strap circumferentially around the lower arm under slight tension. Close the strap by attaching the free attachment module (**3**) to the adjustment module (**4**). Pay attention to comfortable fitting of the Stabilo strap, neither too loose nor too tight. If the Stabilo strap with the easyClick system is not needed, the supporting module (**2**) can be taken from the velour bag.

Removal

After opening the Stabilo strap, take hold of the lower edge of the support and slip the support over the lower arm in the direction of the hand.

Material composition

Polyamide (PA), elastane (EL), silicone (SI)

The sewn in textile label provides the precise material composition.

Cleaning information

 Delicate wash  Do not clean chemically  Do not bleach

 Do not dry in a tumble dryer  Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2019

Medical device |  Single patient – multiple use

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