BORT select®

EpiPlus

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select®

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REF 122 700

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	03
EN	english	Instructions for use	07
FR	français	Mode d'emploi	10
ES	español	Instrucciones de uso	14
IT	italiano	Istruzioni per l'uso	18
NL	nederlands	Gebruiksaanwijzing	22
CS	český	Návod k použití	25
ET	eesti	Kasutusjuhend	28
PL	polski	Instrukcja użytkowania	31
RO	româneșc	Instrucțiunile de utilizare	35

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a circumferential compression elbow compression support with elastic material fabric pads and elastic profile pads.

Indications

Chronic post-traumatic or post-operative soft tissue irritation in the elbow joint area, epicondylitis, tendomyopathy.

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

- In order to avoid damages to the knitted support fabric, only attach the Stabilo strap after successfully determining the size and trial fitting. Goods damaged by the Velcro fastening cannot be returned.

Fitting instructions

Fit the elbow support so that the comfort zone is positioned in the antecubital fossa. Here, observe that you can feel the outer and inner condyle through the hollows in the silicone pad. The circumferential Stabilo tape enables additional individual regulation of the pressure on the silicone pads.

Removal

After opening the circumferential Stabilo strap, take hold of the lower edge of the support and pull it off over the lower arm in the direction of the hand.

Material composition

Polyamide (PA), elastodiene/latex (ED), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash Do not clean chemically Do not bleach

Do not dry in a tumble dryer Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for

the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2019

Medical device |  Single patient – multiple use

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