BORT AchilloStabil Plus Sport

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Sport

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REF 053 100 SP

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	03
EN	english	Instructions for use	07
FR	français	Mode d'emploi	10
ES	español	Instrucciones de uso	14
IT	italiano	Istruzioni per l'uso	17
NL	nederlands	Gebruiksaanwijzing	20
CS	český	Návod k použití	23
ET	eesti	Kasutusjuhend	26
PL	polski	Instrukcja użytkowania	29
RO	româneșc	Instrucțiunile de utilizare	32

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a support for compression of the Achilles tendon with local pressure pads consisting of flexible material. They are intended to have a massaging effect on the Achilles tendon.

Indications

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e.g. in the case of Achilles tendon rupture, bursitis subachillea.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds

- do not use in case of intolerance of one of the materials used
- no re-use this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

- Slip your foot into the donning aid.
- Grasp the upper edge of the support and pull it over your foot.
- Ensure that the profile insert fully covers the Achilles tendon.
- Correct the positioning of the support if applicable.
- When finished, simply pull the slip-on aid forwards out of the support

For additional relief of the Achilles tendon, one of the heel cushions provided can be inserted into the shoe on the side of the ailing leg. In order to avoid leg length differences between the healthy and ailing sides, place the second heel cushion provided in the shoe on the side of the healthy leg.

Removal

To remove the support, please grasp the lower edge and pull it down over the foot.

Material composition

Polyamide/COOLMAX® (PA), polyamide (PA), elastodiene/LYCRA® (ED), silicone (SI).
The sewn in textile label provides the precise material composition.

Cleaning information

 Delicate wash  Do not clean chemically  Do not bleach

 Do not dry in a tumble dryer  Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed

instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2019

Medical device |  Single patient multiple use

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