BORT MalleoStabil-Orthese

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
REF 100 500

D100500I2020-02I003 ML

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DE BORT MalleoStabil-Orthese

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EN BORT MalleoStabil Brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an ankle brace for the stabilisation of the upper ankle joint on one plane. Its function is to restrict pronation and supination via two connected stabilisation elements. The stabilisation elements are attached using adjustable belts.

Indications

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures, ankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, upper ankle joint distortion, fibular ligament rupture, pronation and/or supination trauma.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated

Application risks/important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
 - if wearing at night is prescribed, avoid negative influence on the circulatory system
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use – this medical aid is intended for treating one patient
 - whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

1 When putting the brace on, ensure that the lateral shells are precisely centred with the lower leg. Please set the strap for the sole of the foot to fit the patient's circumstances in order to achieve effective stabilisation. Lay the wide lower strap circumferentially around the shells and attach it with Velcro.

2 Lay the upper strap circumferentially around the shells and attach it with Velcro. If required, adjust the lower strap so that the brace is firmly in position.

3 When the swelling is reduced, it is possible that the MalleoStabil Brace becomes loose. If required, adjust it.

Use of the gel component

4 Remove the gel component from the bottom upwards. When re-attaching the gel component, observe that it is placed directly at its original location. Its correct placing ensures that the plastic edges are covered.

For cold therapy

- Remove the gel component and put it into a deep-freezer for an hour.
- Information:** The gel doesn't harden.
- Re-attach the gel component at its original location in the brace.
- Do not bring the gel component into direct contact with the skin.
- Always wear a sock for best possible comfort.

Foam padding

The gel pads can be exchanged with the foam pads provided as the circumstances require. For this, please be sure to observe the medical prescription.

Use of the foam padding

Remove the foam component from the bottom upwards. When re-attaching the foam pads, observe that they are placed precisely at their original location. Its correct placing ensures that the plastic edges are covered.

Removal

You only need to open the two circumferential belts to remove the BORT MalleoStabil Brace.

Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition

80 % Polypropylene (PP), 20 % polyamide velour (PA)

Cleaning information

- Hand-washing Do not clean chemically Do not bleach
- Do not dry in a tumble dryer Do not iron
- Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fasteners to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the

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