

Summary about Experiment In-Vivo-Gene Editing by CRISPR

Lab workshop and the evaluation of experiments in total done by 20 international student teacher teams guided by mentors from Greece, Czech and Slovakia took us 3 days and showed the following results.

1. Step: Transformation

14 workgroups out of 20 had performed a perfect result. They edited successful the targeted *lac-Z*-gene coding for β -galactosidase on the bacterial chromosome by the active CRISPR/Cas-complex inside the living cells in-vivo. Bacteria colonies of *E. coli* in approach D were existing on the nutrient and their color changed from wild type *lacZ* from blue to white.

6 Workgroups had no colonies in approach D. But other approaches A – C of the 6 workgroups showed expected results (see Fig. 1 and 2).




Fig. 1: Scheme of expected perfect results in the approaches A – D.



Fig. 2: Color of colonies in blue/white screening in the approaches A - D.

2. Step: Confirmation of *LacZ*-Gene Status by Multiplex PCR

PCR amplification of isolated DNA out of blue colonies showed in gel-electrophoresis the length of wild type *lacZ*-fragment of not edited *lacZ*-gene 1100 bp (see Fig. 3).

PCR amplification of isolated DNA out of each investigated white colony of 14 workgroups showed in gel-electrophoresis the length of edited *lacZ*-gene fragment (650 bp).

Each investigated colony showed in gel-electrophoresis the control band for chromosomal DNA with the length of 350 bp (see Fig. 3).




Fig. 3: Band evaluation in gel-electrophoresis after multiplex PCR (M = marker; A, C and D approaches; +c = positive control; -c = negative control).