
STANDARD F C. difficile GDH FIA
STANDARD™ F C. difficile GDH FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer


SPECIMEN COLLECTION AND PREPARATION

- Feces**
- Fecal specimen should be collected in a clean and dry container and collected at any time of the day may be used.
 - Urine should be excluded from the fecal specimen.
 - Fresh fecal specimen should be used in this test. Do not use any transport media to store and transport the fecal specimen.
 - Fecal specimen may be stored in a refrigerator at 2-8°C for up to 24 hours.
 - For prolonged storage, specimens may be frozen under -20°C.
 - It should be brought to room temperature prior to use.

PREPARATION AND TEST PROCEDURE

■ Preparation

- Allow kit components and collected sample to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for STANDARD F C. difficile GDH FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch, and check the test device and the silica gel inside the foil pouch.


If there is no violet colored Check Band on the membrane of the test device, do not use it.

Do not write on the barcode or damage the barcode of the test device.

If the color of silica gel are changed from yellow to green, please do not use the test device.

■ Collection of sample

- Poke the 4 different areas of the fecal sample using the sterile swab to collect a portion of feces (about 40-70mg). In case of liquid sample, soak the sterile swab completely in the liquid fecal sample.

The amount of fecal swab may affect the results. It is required to follow the swab amount of feces as shown in the picture above. Excessive fecal amount may induce a false positive result and slow migration.

- Take off the top of an extraction buffer tube and insert the swab into the tube.

- Swirl the swab at least ten times into the extraction buffer tube to dissolve the sample.
- Remove the swab and dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Tightly press down the filter cap onto the tube.


Do not use the extraction buffer tube without the filter cap. It affects the result value.


■ Analysis of sample

STANDARD TEST™ mode STANDARD F100, F200 and F2400 analyzer

- Prepare a STANDARD F Analyzer and select the 'Standard Test' mode according to the analyzer's manual.
- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.
- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid. Select the sample type on the sample selection screen.
- Apply 3 drops of mixed sample to the sample well in the test device holding the prepared extraction buffer tube upside down.

- After applying the sample, immediately press the 'TEST START' button.
- The analyzer will automatically display the test result within 15 minutes.

- If there is no violet colored Check Band on the membrane of the test device, do not use it.
- Do not put bubbles in the sample well of the test device.
 - Discard the used test kit in proper container, according to your healthcare professional's recommendation.
 - If you feel the test result is inaccurate or do not agree your test result, do not change the treatment and contact your healthcare professional.

INTERPRETATION OF TEST RESULT

| Result | COI (Cutoff index) value | Interpretation |
|--------------------------|----------------------------|--|
| C. difficile GDH-Pos (+) | COI ≥ 1.0 | Positive for C. difficile GDH |
| C. difficile GDH-Neg (-) | COI < 1.0 | Negative for C. difficile GDH |
| Invalid | COI value is not displayed | Retest should be performed with a new test device and a new patient sample |

- Result must be interpreted together with STANDARD F C. difficile Toxin A/B FIA
- False positive result can be caused by intake of antibiotics or other medicine.
- Retest should be performed with newly collected fecal specimen.

EXPLANATION AND SUMMARY

■ Introduction

Clostridium difficile is a gram-positive anaerobe. Infection with *C. difficile* causes severe diarrhea and can be fatal if not diagnosed and treated in a timely manner. *C. difficile* infection is induced in patients by long-term treatment with antibiotics and is commonly found in hospital environment. It is easily transmitted through contact with infected fecal matter. Accurate testing for toxins proved to be difficult due to assays having low sensitivity. Since all strains of *C. difficile* produce large amounts of glutamate dehydrogenase, testing for this antigen has proven to be a better screening tool due to its higher negative predictive value. *Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxicogenic strains, making it an excellent marker for the organism. Diagnostic tests that detect GDH are highly sensitive. They have relatively low specificity for disease-causing *C. difficile* though, because GDH is also produced by non disease-causing *C. difficile* strains and other *Clostridium* species.

■ Intended use

STANDARD F C. difficile GDH FIA is a fluorescence immunoassay for the qualitative detection of *C. difficile* GDH from fecal specimens from persons. This is only for professional use. It provides only an initial screening test result. STANDARD F C. difficile GDH FIA should be used with the appropriate analyzer, STANDARD F Analyzers, manufactured by SD BIOSENSOR.

■ Test principle

STANDARD F C. difficile GDH FIA is based on the immunofluorescence technology with STANDARD F Analyzer, manufactured by SD BIOSENSOR to measure the *C. difficile* GDH in human sample. The specimen from human should be processed for the preparation using the components of STANDARD F C. difficile GDH FIA. After applying the sample mixture to the test device, the complex will be formed on the membrane as the result of the antigen-antibody reaction. The intensity of the fluorescence light generated on the membrane is scanned by STANDARD F Analyzer. STANDARD F Analyzers can analyze the *C. difficile* GDH of the clinical specimen based on a pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- STANDARD F C. difficile GDH FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this instructions before testing.
- STANDARD F C. difficile GDH FIA should be with STANDARD F Analyzer.
- STANDARD F C. difficile GDH FIA should remain in its original sealed pouch until ready to use. Do not use if the pouch is damaged or the seal is broken.
- STANDARD F C. difficile GDH FIA is single use only. Do not re-use it.
- Do not use frozen sample or artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm, soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Check the expiration date printed on the pouch or package.
- Use STANDARD F C. difficile GDH FIA at 15-32°C / 59-60°F.
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

QUALITY CONTROL

■ Calibration
The calibration set test of STANDARD F Analyzers should be conducted according to the analyzers' manual.

[When to use calibration set]

- Before using the analyzer for the first time
- When you drop the analyzer
- Whenever you do not agree with your result
- When you want to check the performance of an analyzer and test device

[How to use calibration set]

- Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.
- Select the 'Calibration' menu.
 - The specific calibration set is included with the analyzer.
 - Insert the CAL-1 first, and then insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.
 - STANDARD F analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in 'Standard Test' mode. If 'EEE' message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the 'EEE' message still appears.

■ External quality control

Quality control testing should be run to check the performance of STANDARD F C. difficile GDH FIA and STANDARD F Analyzers. STANDARD F C. difficile GDH Control manufactured by SD BIOSENSOR should be used for quality control testing. Control test should be conducted in accordance with the instructions of STANDARD F C. difficile GDH Control.

- Control test should be run:
- once for each new lot.
 - once for each untrained operator.
 - as required by test procedures in instructions for use of STANDARD F C. difficile GDH Control and in accordance with local, state and federal regulations or accreditation requirements.

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PERFORMANCE CHARACTERISTICS

■ Clinical evaluation

- Limit of Detection
The limit of detection of STANDARD F C. difficile GDH FIA is as follows;

| Source | Concentration |
|--------------|---------------|
| Bionote-U156 | 1.34 ng/ml |

- Precision
STANDARD F C. difficile GDH FIA totally meet the acceptance criteria of repeatability and reproducibility.

| Panel No. | Operator | Operator 1 | | | Operator 2 | | | Operator 3 | | |
|-----------|-------------------|------------|-------|-------|------------|-------|-------|------------|-------|-------|
| | | Site 1 | | | Site 2 | | | Site 3 | | |
| | | Lot | Lot 1 | Lot 3 | Lot 1 | Lot 2 | Lot 3 | Lot 1 | Lot 2 | Lot 3 |
| C | CIFF-GDH-1801-P02 | 100% | 100% | 100% | 100% | 100% | 100% | 100% | 100% | 100% |
| | CIFF-GDH-1801-P01 | 97% | 95% | 97% | 97% | 98% | 98% | 97% | 95% | 97% |
| | CIFF-GDH-1801-N00 | 0% | 0% | 0% | 0% | 0% | 0% | 0% | 0% | 0% |

- Relative sensitivity & specificity (LFA product compared)
Relative sensitivity: 97.5% (80/82) [95% CI: 87.23 - 100%]
Relative specificity: 100% (77/77) [95% CI: 97.72 - 100%]
- Cross Reactivity Study
STANDARD F C. difficile GDH FIA did not show any sign of cross reactivity with the concentrations list of each cross-reacting microorganisms.

| Potential Cross-Reactant | Concentration |
|--------------------------|----------------------------|
| Adenovirus type 1 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 2 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 3 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 4 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 5 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 7 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 11 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 18 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 23 | 2.0x10 ⁶ cfu/mL |

| | |
|-------------------------------|----------------------------|
| Adenovirus type 40 | 2.0x10 ⁶ cfu/mL |
| Adenovirus type 55 | 2.0x10 ⁶ cfu/mL |
| Bacteriodes fragilis | 2.0x10 ⁶ cfu/mL |
| Candida albicans | 2.0x10 ⁶ cfu/mL |
| Escherichia coli | 2.0x10 ⁶ cfu/mL |
| Klebsiella pneumoniae | 2.0x10 ⁶ cfu/mL |
| Peptostreptococcus anaerobius | 2.0x10 ⁶ cfu/mL |
| Pseudomonas aeruginosa | 2.0x10 ⁶ cfu/mL |
| Staphylococcus aureus | 2.0x10 ⁶ cfu/mL |
| Staphylococcus epidermidis | 2.0x10 ⁶ cfu/mL |

- Interference Substances
STANDARD F C. difficile GDH FIA is not affected by below substances.

| Potential Interfering Substance | Concentration |
|---------------------------------|---------------|
| Barium sulfate | 5% w/v |
| Human whole blood | 40% v/v |
| Hemoglobin | 12.5% w/v |
| Imodium | 5% w/v |
| Pepto-Bismol | 5% w/v |
| Vancomycin | 0.2% w/v |

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BIBLIOGRAPHY

- Wren, M.W.D, et al. "Laboratory diagnosis of infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
- Vaishnavi, Ch., "Clinical spectrum & pathogenesis of difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499.
- Sharp SE, Ruden LD, Pohj J, et al. Evaluation of the C.diff Quik Chek Complete Assay, a new glutamate dehydrogenase and A/B toxin combination lateral flow assay for use in rapid, simple diagnosis of Clostridium difficile disease. J Clin Microbiol. 2010;48:2082-2086.
- Surawicz CM, Brandt LJ, Binion DG, et al. Guidelines for diagnosis, treatment, and prevention of Clostridium difficile infections. Am J Gastroenterol. 2013;108:478-498.
- Reiler ME, Akcabas RC, Lema CA, et al. Comparison of two rapid assays for Clostridium difficile common antigen and a C difficile toxin A/B assay with the cell culture neutralization assay. Am J Clin Pathol. 2010;133:107-109.

SYMBOL

| REF | Reference number | Caution | Use by | LOT | Batch code | Consult Instructions for Use | Do not re-use |
|-----|----------------------|---------|--------------|-----|---------------------|---|-------------------------|
| IVD | In vitro Diagnostics | Note | Manufacturer | | Date of manufacture | Contains Sufficient for ∞ Tests | Keep away from sunlight |
| | | | | | | | |