EN

STANDARD F CK-MB FIA

STANDARD™ F CK-MB FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF

F-CKB

EXPLANATION AND SUMMARY

[Introduction]

Creatine Kinase (CK), also known as Creatine Phosphokinase or Phospho-creatine Kinase is an enzyme expressed by various tissues and cell types. Distribution of cell membrane CK isoenzyme is CK from the cardiac cytosolic membrane. Creatine kinase is a dimer formed by the association of two subunits from muscle tissue type or nerve cells (brain type). These subunits associate to form three isoenzyme forms: CK-BB, CK-MM and CK-MB. These isoenzymes are expressed at different levels in various human tissues. Though CK-MM is the most abundant CK isoenzyme in the cardiac muscles, CK-MB constitutes about 20% of the total CK in the healthy heart muscle tissue. Elevated level of total CK may indicate myocardial damage, but it can be observed in patients with skeletal muscle injury or certain other disorders but as CK-MB is more specific to myocardial tissue, CK-MB levels along with total CK can be considered as an important diagnostic indicator of myocardial infarction. The concentration of CK-MB in the healthy adult myocardium is low concentration but it shows greater increase in case of myocardial infarction, myocarditis, myocardial edema and infarction. CK-MB is known to be more sensitive and earlier marker of myocardial infarction because it has a lower basal level and a narrower normal range. Medical literature commonly reveals that following an acute myocardial infarction, CK-MB levels become elevated in 4 to 9 hours after the onset of chest pain, attain peak levels in 24 to 48 hours and return to baseline within 2 to 3 days. Thus, measurement of CK-MB in serum or whole blood is the most important clinical application of CK measurements in clinical chemistry. STANDARD F CK-MB FIA, employing immunofluorescent detection system with STANDARD F Analyzer, provides significantly fast, easy and accurate system to identify CK-MB in serum and whole blood as an aid in the diagnosis of myocardial infarction and enables supportive treatment decisions.

[Intended use]
STANDARD F CK-MB FIA is a fluorescence immunoassay for the quantitative determination of total Creatine Isoenzyme MB (CK-MB) levels in human serum and whole blood using STANDARD F Analyzers, manufactured by SD BIOSENSOR. This test is an in vitro diagnostic use and intended for use as an aid in the screening and monitoring of myocardial infarction (MI).

[Test principle]
When a human specimen is applied into the sample well of the test device, the specimen migrates through the membrane from the sample well. If Creative Kinase Isoenzyme-MB (CK-MB) is present, it will be bound to the detector coupled to a europium microparticle through its antibody binding site. The europium microparticle complexed with the capture antibody on the test line where it is detected by STANDARD F Analyzer. If CK-MB is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the CK-MB in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

[Kit contents]

① Test device ② Fixed volume dropper (100 μl) ③ Disposable dropper (100 μl) ④ Extraction buffer

⑤ Instructions for use

• STANDARD F Analyzer

KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F, out of direct sunlight. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

1. STANDARD F CK-MB FIA is for in vitro diagnostics use only.
2. Carefully follow all instructions described in these instructions before testing.
3. STANDARD F CK-MB FIA should be used with STANDARD F Analyzer.
4. STANDARD F CK-MB FIA should remain in its original sealed pouch until ready to use. Do not use if the pouch is damaged or the seal is broken.
5. STANDARD F CK-MB FIA is single use only. Do not re-use.
6. Do not use frozen samples.
7. Do not use any artificial materials.
8. Place the analyzer on a flat surface when in use. Wash your hands in warm, soapy water. Rinse well and dry completely before testing.
9. Do not use any reagent not according to the proper method.
10. Silica gel in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating silica gel beads change from yellow to green, the test device in the pouch should be discarded.
11. Use a fixed volume dropper (100 μl) - [3] labeled on the pouch for sample extraction purpose only. Do not use a disposable dropper.
12. Use a disposable dropper (100 μl) - [2] labeled on the pouch for mixing the sample, and for dispensing the sample mixture into the test device.
13. Check the expiration date on the pouch or package.
14. Check the volume (100 μl) of extraction buffer.
15. Use the STANDARD F CK-MB FIA at room temperature.
16. All kit components must be at room temperature 30 minutes before running the assay.
17. Do not write on the bar code or damage the bar code of the test device.

SPECIMEN COLLECTION AND PREPARATION

[Whole blood]

1. Collect the whole blood into the commercially available EDTA tube by venipuncture.
2. It is recommended that collected venous whole blood samples are used immediately. If venous whole blood in an anticoagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 8 hours after collection.
3. Do not use hemolyzed blood samples.

[Serum]


1. Collect the whole blood into the commercially available plain tube, NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
2. Serum samples may be stored at room temperature for up to 1 day or at 2-8°C/36-46°F for up to 3 days prior to testing.
3. For over 3 days storage, specimens may be frozen under -20°C/-4°F for up to 1 month.
4. It should be brought to room temperature prior to use.

TEST PROCEDURE

[Preparation]

1. Allow kit components and collected sample to room temperature at least 30 minutes before starting the test.
2. Carefully read instructions for the STANDARD F CK-MB FIA.
3. Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.

4. Open the foil pouch, and check the test device and the silica gel pack inside the foil pouch.


CAUTION: Do not write on the bar code or damage the bar code of the test device.

If the color of moisture indicators are changed from yellow to green, please do not use the test device.

[Analysis of sample]

1. Using a 'STANDARD TEST' mode
2. STANDARD F Analyzer and select the 'Standard Test' mode according to the analyzer's manual.
3. In case of STANDARD F200 and F2400 analyzer, go to the 'Workplace' in the main screen. And select the 'Run Test'.
4. Take the test device out of the foil pouch.

9. The analyzer will automatically display the test result within 10 minutes.

OK START Test Start <F100> <F200> <F2400>

10. Read the test result after 10 minutes.

• Do not put bubbles in the sample well of the test device.
• Discard the used test kit in proper container, according to your healthcare professional's recommendation.
• If you feel the test result is inaccurate or do not agree your test result, do not change the treatment and contact your healthcare professional.

INTERPRET THE TEST RESULT

- Results should be considered in conjunction with the clinical history and other data available to the physician.
- If an error message appears on the analyzer's screen, refer to the analyzer's manual.

Measuring Range 1-200 ng/mL

Out of measuring range - If below 1 ng/mL, '1' message will be displayed.
- If above 200 ng/mL, '200' message will be displayed.

Cut off value for MI <5.0 ng/mL

* People with skeletal muscle damage may have elevated CK-MB levels. To help differentiate heart attack from skeletal muscle damage, a CK index can be calculated using CK-MB and total CK as follows: CK index = (CK-MB, ng/mL) x 100 / total CK activity, IU/L. If CK-MB is elevated and the CK index is higher than 2.5 to 3.0, heart damage is likely.

NOTE: The CK-MB reference ranges are provided for orientational purpose only. Clinicians should use the test results in conjunction with the patient's other diagnostic findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

QUALITY CONTROL

The calibration set of STANDARD F Analyzers should be conducted according to the analyzers' manual.

[When to use calibration set]

1. Before using the analyzer for the first time
2. When you drop the analyzer
3. Whenever you receive with your result
4. Whenever you check the performance of an analyzer and test device

[How to use calibration set]

Calibration set is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

• The specific calibration set is included with the analyzer.

3. Insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.

- The STANDARD F analyzer automatically calibrates and identifies the optical performance through measuring the membrane of the test device whenever the test is conducted in Standard Test mode. If 'EEE' message displays on the screen, means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the 'EEE' message still appears.

EXTERNAL QUALITY CONTROL

Quality control testing should be run to check the performance of STANDARD F CK-MB FIA and STANDARD F Analyzers. STANDARD F CK-MB Control manufactured by SD BIOSENSOR should be used for quality control testing. The test should be conducted in accordance with the instructions of STANDARD F CK-MB Control.

• Once for each new lot.

- once for each untrained operator.

- as required by test procedures in instructions for use of STANDARD F CK-MB Control and in accordance with local, state and federal regulations or accreditation requirements.

BIBLIOGRAPHY

1. Mier J, Andrade-Orozco E, Dienst A, et al. Early detection of acute myocardial infarction by measurement of mass concentration of creatine kinase-MB. Am J Cardiol 1991;68:1545-1550.
2. Boucek RJ Jr; Kasselberg AG; Boerth RC; Parrish MD; Graham TP Jr (July 1982). "Myocardial injury in infants with congenital heart disease: evaluation by creatine kinase MB isoenzyme analysis". Am J Cardiol. 50 (1): 129-132. doi:10.1016/0002-9140(82)90705-7. PMID 7102020.
3. "Use of Cardiac Markers in the Emergency Department". eMedicine. Retrieved 2009-01-10.
4. Guyz PM (December 1977). "Creatine phosphokinase-MB (CPK-MB) and the diagnosis of myocardial infarction". West J Med. 127 (6): 455-60. PMC 1237914 Freely accessible. PMID 339548.

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