EN

REF FPNE01G

S. pneumoniae Ag FIA

STANDARD™ F S.pneumoniae Ag FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

STANDARD™

EXPLANATION AND SUMMARY

[Introduction]

S. pneumoniae (*S. pneumoniae*) is a gram-positive and catalase-negative organism, also called pneumococcus since it causes the pneumococcal infection. *S. pneumoniae* is involved in upper respiratory infections such as otitis media and sinusitis as well as diseases of the lungs as pneumonia, bacteremia and meningitis. *S. pneumoniae* is spread through indirect or direct contact with an *S. pneumoniae* carrier and spread of almost all the upper respiratory infections are due to the carriage by *S. pneumoniae* bacteria. Especially, rapid and reliable diagnosis of *S. pneumoniae* is an important part of the control process (CAP) of community carriers. The collection of samples for analysis and proper transport (CAP) is essential to reduce the risk of pneumonia. Bacterial meningitis is serious and it can be deadly. Most people recover from meningitis. However, permanent disabilities such as hearing damage, hearing loss and learning disabilities can result from infection. Antibiotic treatment for bacterial meningitis, it is important to start treatment as soon as possible.

[Intended use]

STANDARD F *S. pneumoniae* Ag FIA is a fluorescence immunoassay for the qualitative detection of *S. pneumoniae* CPS Ag in urine samples. The sample is taken with a pipette with respiratory symptoms and in cervical spinal fluid (CSF) samples from patients with meningitis. The test is used for professional diagnosis and intended as aid to early diagnosis of *S. pneumoniae* infection. It provides only an initial screening test result. Specific alternative diagnosis method should be performed in order to obtain the confirmation of *S. pneumoniae* infection.

[Test principle]

STANDARD F *S. pneumoniae* Ag FIA is based on immunofluorescence technology with STANDARD Analyzer to detect *S. pneumoniae* antigen. STANDARD F *S. pneumoniae* Ag FIA has control line and test line which is coated with monoclonal anti-PPR and rabbit anti-S-pneumoniae IgG each. The patient's sample is applied into the sample well of the test device and the sample migrates through the membrane. If *S. pneumoniae* antigen is present, it reacts with enzyme conjugated anti-S-pneumoniae IgG and anti-S-pneumoniae CPS in the conjugation pad and form antibody-antigen fluorescence particle complexes. These complexes move along the membrane to be captured by the rabbit anti-S-pneumoniae IgG antibodies on the test line. The intensity of the fluorescence of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer manufactured by SD BIOSENSOR. STANDARD F Analyzer can analyze the presence of the analyte in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

[Kit contents]

① Test device ② Fixed volume pipette (100μl) ③ Positive control ④ Negative control

⑤ Control extraction buffer ⑥ Instructions for use

[Materials required but not provided]

• Timer

KIT STORAGE AND STABILITY

Store the kit at 2-30°C/56-86°F, out of direct sunlight. Kit materials are stable until expiration date printed on outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

1. Do not eat or drink.
2. Use the STANDARD F *S. pneumoniae* Ag FIA at 15-32°C / 59-90°F and 10-90%RH.
3. Do not use the kit if the pouch is damaged or the seal is broken.
4. Do not smoke, drink or eat while handling specimens.
5. Wear appropriate personal equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after you perform.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Handle all specimens as if they contain infectious agents.
8. Observe standard precautions against microbiological hazards throughout testing procedures.
9. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory cleanup and disposal must be handled and discarded in accordance with local, state, and national regulations.
10. Silica gel in the pouch is an absorbent moisture and keep humidity from affecting products. If the moisture indicating silica gel beads change from yellow to green, the test device in the pouch should be discarded.
11. The test device of the kit is used by analyzer to identify the type of test being run and to identify the individual test device so as to prevent a second read of the test device by another operator.
12. After use, do not use the test device after taking out of a foil pouch.
13. As the detection reagent is a fluorescent compound, no visible results will form on the test device. The STANDARD F Analyzers authorized by SD BIOSENSOR must be used for result interpretation.
14. Improper specimen collection, handling or transport may yield inaccurate results.
15. Do not write on the bar code or damage the bar code of the test device.
16. If the test result with positive/negative control is abnormal, do not use the kit.

SPECIMEN COLLECTION AND PREPARATION

[Urine]

1. Urine specimens should be collected in a clean and dry container and collected at any time of the day may be used.
2. Urine specimen may be stored at room temperature for up to 24 hours or at 2-8°C / 36-46°F for up to 1 week prior to testing.
3. For over 1 week storage, specimen may be frozen and stored at -40°C/-40°F.
4. Frozen specimens should be thawed and fully mixed by vortexing before testing.
5. If the specimen exhibits visible precipitates, it should be allowed to settle to obtain a clear supernatant for testing.

[CSF]

1. CSF specimen should be collected in a dry and sterilized standard container and collected at any time of the day may be used.
2. Specimen may be stored at room temperature for up to 24 hours prior to testing.
3. When the specimen stored at 2-8°C condition, the test should be performed within anytime up to 1 week.
4. If the specimen contains visible precipitates, it should be allowed to settle to obtain a clear supernatant for testing.

TEST PROCEDURE

[Preparation]

1. Allow kit components and collected urine sample to room temperature (15-30°C / 59-86°F) a minimum of 30 minutes before use.
2. Carefully read instructions for using the STANDARD F *S. pneumoniae* Ag FIA.
3. Check the expiry date at the back of the foil pouch. Use another lot, if expiry date has passed.

4. Open the foil pouch, and check the test device and the silica gel pack in the foil pouch.

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
- [BIBLIOGRAFÍA]**
1. Streptococcus pneumoniae colonization: the key to pneumococcal disease, THE LANCET Infectious Diseases, Vol. 4, No. 3, p144-154, March 2004.
 2. Antibody to Antigenic Polysaccharides of Streptococcus pneumoniae - Prevalence, Persistence, and Response to Recuperation, Clin Infect Dis, 1993; 17(1):66-73.
 3. Emergence of Antibody to Capsular Polysaccharides of Streptococcus pneumoniae During Outbreaks of Pneumonia, Clin Infect Dis, 1993; 17(1):66-73.
 4. Spread of Drug-Resistant Streptococcus pneumoniae in Asian Countries: Asian Network for Surveillance of Resistant Pathogens (ANSOP) Study, Clinical Infectious Diseases, 1999; 28:1206-11.

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