

# Analysis of difficult samples with ICP-OES

Agilent 5110 ICP-OES

Sebastien Sannac  
*Atomic Spectroscopy Product Specialist*



# Benefits of the 5100



## Lowest cost of ownership

- Fastest sample throughput
- Low gas consumption

## Enhanced Performance

- Analytical performance
- System robustness and reliability

## Simple Operation

- Hardware
- Software

# Difficult Samples



# Difficult Samples

Running highly concentrated, complex samples brings a range of new challenges:

## Solid Build Up in Torch Injector

- Changes gas flow velocity and observation zone. This manifests as drift.

## Solid build up in Nebulizer

- Blockage leads to poor precision and signal drift.

## Organic solvent based samples


- Precision can be impacted.

## Interferences from the sample matrix on the analytes


- Results accuracy and truness can be impacted

# Difficult Samples

Samples with high levels of dissolved solids can start to deposit in the injector.




Vertical Torch – 25% TDS capability



Horizontal Torch – 5% TDS capability


# Difficult Samples

5110 ICP-OES series : Vertical Torch System and Dual View capacity



- Vertical torch for matrix tolerance.
- Axial view for high sensitivity.
- Radial view for robust measurement (EIEs).
- DSC technology allows axial and radial light to be detected simultaneously.

# Difficult Samples



- Desolvation of high TDS sample occurs due to the fast flowing gas through the small nebulizer orifice – leads to crystalline material depositing
- High Total Dissolved Solids samples can lead to poor precision.

White salt crystal lodged in nebulizer orifice





# Difficult Samples

## *Argon Humidifier Accessory (AHA)*



- Permeable PTFE membrane allows H<sub>2</sub>O vapor into the nebulizer gas.
- Humidifies nebulizer gas to help reduce:
  - Nebulizer blockages owing to salt build-up.
  - Long term **drift** from deposits in torch injector.





# Difficult Samples


## *Double-Pass Spraychamber*

- Double-pass spraychambers reduce sample loading on plasma:
  - Effectively sorts out larger droplets that destabilize plasma and degrade precision.
  - Reduces build-up of deposits in torch injector that can change nebulizer gas velocity through plasma central channel causing long term drift.



# 250 ppb Multi-element in 25% NaCl

*<2.4% RSD Long Term Stability with VDV Configuration (axial only mode)*



OneNeb, Double pass glass cyclonic, Blk/Blk uptake pump tubing, Ar humidifier, 1.8mm injector demountable torch, 15rpm peristaltic pump

# Organic Samples

- Select appropriate sample introduction system.
- Incorrect selection leads to plasma instability – poor precision and even plasma extinguishing.
- Typically requires a narrower Internal Diameter injector than standard e.g.
  - 1.4 mm ID torch injector for non-volatile organics such as kerosene.
  - 0.8 mm ID torch injector for volatile organics such as gasoline or naphtha.
- Use appropriate solvent resistant pump tubing.



# Organic Samples



- For volatile organics use a cooled spraychamber (see below) at sub-ambient conditions e.g.  $-10^{\circ}\text{C}$
- Organic resistant pump tubing such as marprene, solvent flex and viton.



May need to use some oxygen addition into auxillary Ar flow to prevent carbon deposition on injector




# Fitted Background Correction (FBC)




- Fully automated background correction.
- Automatically models the background signal under the analyte peak using a sophisticated mathematical algorithm.
- Traditional Off-Peak background correction requires the background points to be manually selected.

# Fitted Background Correction (FBC)




- FBC calculates the true background signal – improving accuracy.
- Off Peak Left and Right background correction needs to be adjusted manually for the sample shown on the bottom left.
- FBC requires no user adjustment unlike Off Peak background correction.



# Removing Spectral Interferences

## FACT – Fast Automated Curve- Fitting Technique

Possible interferences on Cd (226.502 nm)



- Resolves extremely complex spectral interferences.
- Allows resolution of interferences as close as 3 pm.
- Re-analysis of samples not required.
- Quick and Easy to set up.



# Removing Spectral Interferences


## FACT – Fast Automated Curve- Fitting Technique

Possible interferences on Cd (226.502 nm)

Symbol	Wavelength (nm)	Ion	Intensity
Ni	226.446	II	3752.8
Nb	226.455	II	122.8
Fe	226.459	II	16.0
Ir	226.461	I	98.0
Os	226.465	I	10.4
Ta	226.465	II	1.6
Ru	226.470	I	3.3
Mo	226.473	II	4.8
Co	226.488	I	10.6
Mo	226.495	II	2.5
<b>Cd</b>	<b>226.502</b>	<b>II</b>	<b>38378.5</b>
Fe	226.505	I	4.7
Ti	226.510	II	7.8
Mo	226.510	II	1.5
Ir	226.516	II	112.1
Cr	226.519	II	0.6
Nb	226.521	II	0.9
Co	226.526	II	2.1
W	226.533	II	78.4
Ni	226.535	II	12.5

- Resolves extremely complex spectral interferences.
- Allows resolution of interferences as close as 3 pm.
- Re-analysis of samples not required.
- Quick and Easy to set up.

# FACT – How Does it Work?



- Uses a peak modelling approach.
- Models the background, analyte and the interference.
- Deconvolutes the analyte peak from nearby interference peaks.
- When a model is created it can be stored in the library and used again.

# FACT – In Practice

## FACT Models for Interference Correction

Number of interferents:

Edit model names  ⓘ

Show library  ⓘ

ⓘ

### Worksheet models

	Label (Wavelength nm)	Blank	Matrix	Analyte	Interferent 1	Test
▶	Cd (226.502)	Deionised Water Blank	STD 0 (Blank)	STD 1 (0.1 ppm w/v)	500 ppm w/v Fe	

- Select the solution to use as a blank.
- Select the solution to use as a matrix.
- Select the solution to use as the analyte.
- Select the solution to use as the interferent.

# FACT – In Practice

## Without FACT

Solution Label	Cd 214.439 nm ppm	Cd 226.502 nm ppm
1-Rcd 1 280514-1/1	0.095	0.122
2-14E029905-021/1	0.030	0.055
3-14E029872-011/1	0.010	0.047
4-14E029905-002/1	0.042	0.065
5-14E029790-002/1	0.009	0.030
6-14E029896-002/1	0.012	0.035
7-14E029102-006/1	0.010	0.033
8-14E029896-011/1	0.019	0.087
9-14E029488-002/1	0.011	0.035
10-14E029943-001/1	0.010	0.057
11-14E029896-012/1	0.027	0.105
12-14E030194-025/1	0.004	0.024
13-14E029605-005/1	0.026	0.119

- FACT model saved to the library for future use.
- IEC can still be used, however, this requires the interferences to be known before the analysis.
- In conjunction with IntelliQuant this can be a very powerful and quick spectral deconvolution tool.

# FACT – In Practice

## With FACT

Solution Label	Cd 214.439 nm ppm	Cd 226.502 nm ppm
1-Rcd 1 280514-1/1	0.089	0.089
2-14E029905-021/1	0.024	0.024
3-14E029872-011/1	0.001	0.001
4-14E029905-002/1	0.036	0.036
5-14E029790-002/1	0.004	0.004
6-14E029896-002/1	0.006	0.005
7-14E029102-006/1	0.004	0.004
8-14E029896-011/1	0.004	0.005
9-14E029488-002/1	0.005	0.004
10-14E029943-001/1	0.000	0.001
11-14E029896-012/1	0.008	0.008
12-14E030194-025/1	0.000	0.001
13-14E029605-005/1	0.005	0.004

- FACT model saved to the library for future use.
- IEC can still be used, however, this requires the interferences to be known before the analysis.
- In conjunction with IntelliQuant this can be a very powerful and quick spectral deconvolution tool.



# Summary

Vertical torch reduces blockages.  
Argon humidifier to avoid nebulizer clogging.



Dedicated sample introduction system (S/C, torch),  
Control of the S/C temperature



Fitted and FACT background correction ensure the  
right results are reported.



IntelliQuant gives extra sample information –  
excellent for determining interferences.



# The Agilent Atomic Spectroscopy Portfolio



Agilent's **55 and 200 Series** includes *Fast Sequential flame AA* and *high performance furnace*.



Agilent's **4210 MP-AES** runs on air for the *lowest cost of ownership* and *improved safety*.



Agilent's **5110 ICP-OES** includes the *world's most productive*, and *only Synchronous Vertical Dual View ICP-OES*.



Agilent's **7800 & 7900 ICP-MS** are *robust, sensitive, accurate, and easy to use quadrupole ICP-MS*.




Agilent's **8900 ICP-MS** with **MS/MS mode** provides *unique control of interference removal in reaction mode*.

Full portfolio to cover your elemental analysis needs



Agilent Technologies

Tack!



**QUESTIONS?**

*[sebastien\\_sannac@agilent.com](mailto:sebastien_sannac@agilent.com)*



Agilent Technologies