

Routine Trace Metal Analysis of Marine Fuels Using Microwave Acid Digestion

Multiwave PRO

Dr. Markus Michaelis
markus.michaelis@anton-paar.com




Anton Paar – Company Profile

Establishment	1922
Owner	Charitable Santner Foundation
CEO	Dr. Friedrich Santner
Chairman of the Board	DI Ulrich Santner
Number of employees (*as of the end of 2015)	967* Headquarters 2.259* worldwide
Turnover 2015	~ 260 million Euros
Export share	93 %
Investment in R&D	20 % of the annual turnover



Company Presence - Worldwide



- ▶ .. to define the quality
- ▶ .. to prevent deactivation of catalysts
- ▶ .. to monitor and analyze contaminations
- ▶ ...to follow regulations and much more...

Common sample preparation techniques:

- ▶ Dry ashing with open vessel acid digestion (e.g. IP 501, ASTM D5708-B, or D5863-A)
- ▶ Direct dilution using organic solvents (e.g. ASTM D4951, D5708-A, D5863-B, D5185)
- ▶ Microwave-assisted acid digestion (ASTM D7876)



Trace metal determination
ICP-OES, ICP-MS, FAAS

- ▶ Dilution Method:
 - ▶ Metal particles $> 5 \mu\text{m}$ cannot be measured with the dilution method
 - ▶ Clogged nebulizers
 - ▶ Employed organic solvents interfere with modern measuring systems
 - ▶ Special measuring equipment required

- ▶ Dry Ashing & Open vessel digestion
 - ▶ Losses of volatile elements
 - ▶ Increased risk of contamination
 - ▶ Time consuming (in the range of 6 hrs to 1.5 days)

► Multiwave PRO - The Master of Sample Preparation

► Possible Applications:

- **Acid Digestion / Leaching**
- Oxygen Combustion
- UV – Digestion
- Protein Hydrolysis
- Solvent Extraction
- Acid Evaporation
- Drying
- Synthesis



- ▶ Shortest process times
- ▶ Reliable acid digestions up to 300 °C @ 80 bar (1160 psi)
- ▶ No losses of volatiles & minimized risk of contamination
- ▶ Minimized reagent consumption
- ▶ Full reaction control:
 - ▶ Wireless temperature control in each position
 - ▶ Wireless pressure control in each vessel
- ▶ Easy data transfer via USB, LAN
- ▶ Safety – CE, GS & NRTL approval



Why High Temperature? – Complete Digestions

Samples: 300 mg Lubrication Oil



170 °C

190 °C

230 °C


260 °C

Same holding time, different temperatures

**Higher temperature
and pressure →
Better digestion
result!**

High Performance Rotor 8N

- ▶ 8 vessels for high performance
 - ▶ Unique pressure sensor system for simultaneous control of all reaction vessels
 - ▶ Wireless data transmission
 - ▶ Robust rotor upper plate
 - ▶ Temperature control on all positions (via IR sensor)
 - ▶ T-sensor in one reference vessel (optional)
-
- ▶ **XF100: 260 °C @ 60 bar**
 - ▶ **XQ80: 300 °C @ 80 bar**
 - ▶ p_{MAX} : 140 bar




Quick & simultaneous

- ▶ Measures pressure (p)
 - ▶ Measures pressure increase rate
- } → Every 20 milliseconds


If limits are exceeded:

- ▶ MW Power is reduced
- ▶ Cooling is increased



High-speed pressure sensing in action

XQ Digestion with spontaneous reactions



- ▶ @ 10 + 12 min: p-jump of 35 bar
- ▶ Controlled by power cut and cooling up
- ▶ No indication of vessel number, probably different vessels
- ▶ Normal process continuation afterwards

▶ ASTM D7876

- ▶ Standard Practice for Sample Decomposition Using Microwave Heating (With or Without Prior Ashing) for Atomic Spectroscopic Elemental Determination in Petroleum Products and Lubricants.

Applicable to:

- ▶ Petroleum products and lubricants
 - ▶ Greases
 - ▶ Additives
 - ▶ Lubricating oils
 - ▶ Gasolines
 - ▶ Diesels
- ▶ Fossil fuel products:
 - ▶ Coal
 - ▶ Fly ash
 - ▶ Coal ash
 - ▶ Coke
 - ▶ Oil shale

- ▶ ASTM D4951, ASTM D5185, ASTM D7876

- ▶ 3 Samples:

Engine Oil Lubricant, Crude Oil, Residual Fuel Oil



ASTM Interlaboratory Crosscheck NIST SRM 1634c

- ▶ Multiwave PRO equipped with Rotor 8NXF100


- ▶ Two digestion reagents:

- 1.) 7 mL HNO₃ (65%) + 1 mL H₂O₂ (30%)

- 2.) 7 mL HNO₃ (65%) + 1 mL H₂O₂ (30%) + 200 μL of a 50% NH₄F

- ▶ Solutions measured with ICP-OES (1+2) and with ICP-MS (2)

Digestion Program of Crude Oil Samples



Vessel surface temperatures
~ 200 °C

Digestion conditions:
 $T > 200\text{ °C}$ & $p \sim 60\text{ bar}$

Power is reduced as p-Limit
is reached

Pressure of 60 bar

Residual Fuel Oil (NIST SRM 1634c)

	D7876 ICP-OES [$\mu\text{g}\cdot\text{g}^{-1}$]	D7876 ICP-MS [$\mu\text{g}\cdot\text{g}^{-1}$]	Certified value [$\mu\text{g}\cdot\text{g}^{-1}$]
Al	3.6 ± 0.6		
Co		0.14 ± 0.05	0.1510 ± 0.0051
Fe	49.5 ± 1.5	47.8 ± 2.9	
Mg	2.0 ± 0.1		
Na	42 ± 2	33 ± 2	(37)*
Ni	15.7 ± 0.2	17.2 ± 0.3	17.54 ± 0.21
V	27.1 ± 0.4	28.6 ± 0.2	28.19 ± 0.40

ICP-OES: n = 6 (mean of digestions with and without NH_4F)

ICP-MS: n = 3 (only digestions with NH_4F)

*...information value


	2-step, 0.6 g [$\mu\text{g}\cdot\text{g}^{-1}$]	2-step + NH_4F 0.6 g [$\mu\text{g}\cdot\text{g}^{-1}$]	Referencevalue [$\mu\text{g}\cdot\text{g}^{-1}$]
Al	8.1 ± 0.6	8.1 ± 0.6	7.5 ± 1.1
Si	10.2 ± 0.3	4.2 ± 2.2	15 ± 2.2
Ca	48.7 ± 27.6	23.8 ± 0.9	13 ± 2.5
Fe	86.1 ± 0.8	76.7 ± 7.0	53 ± 12
Ni	39.3 ± 6.7	37.2 ± 6.7	39 ± 2.7
Na	26.8 ± 1.4	25.8 ± 0.9	22 ± 4.3
V	80.8 ± 9.3	78.5 ± 5.6	78 ± 4.3

- ▶ Good agreement for most elements
- ▶ Si: problem with sampling heterogeneity?
- ▶ Ca, Fe contamination issues

- ▶ Measured values are in a good agreement
- ▶ Values from ICP-OES and ICP-MS are well comparable → reliable sample preparation procedure
- ▶ Low standard deviations for the 6 samples → 0.5 g of sample are sufficient to achieve representative results.
- ▶ Digestion time reduction from typically **6-8 hours to < 90 minutes**
- ▶ Time for 2-step procedure less than 120 min incl. handling



Pre-reaction

1st Digestion: 0.6 g Fuel Oil 463



2nd step

2nd Digestion of 0.6 g Fuel Oil 46



- ▶ Highest temperatures (300 °C) → lowest interferences during measurement
- ▶ Complete digestions of demanding samples
- ▶ Closed vessel digestion: No losses of volatiles & minimized risk of contamination
- ▶ Fast sample preparation method
- ▶ Minimized reagent consumption
- ▶ Full reaction control for highly reactive samples

