

# VIROMER® ONE RED

## Quick Guide

plasmid and mRNA transfection



ONE step  
transfection

### Recommendations

**Viromer® ONE RED** is fully compatible with cell culture media, sera or antibiotics. Seed cells in complete medium the day before transfection. For good practice, replace medium before transfection.

**Cell culture and plating:** grow cells to reach 60-80% of confluency at the time of transfection.

Use volumes and cell numbers as per the table below. Optimal density for transfection is highly variable within cell types and need to be determined empirically.

Multiwell plate type	96	24
<b>Adherent cells</b>		
Cells seeded per well	12,000	60,000
*Range	± 3,000	± 20,000
<b>Suspension cells</b>		
Cells seeded per well	48,000	240,000
*Range	± 12,000	± 80,000
Medium per well	0.1 ml	0.5 ml

\* in reverse transfection protocols, cell numbers should be on the higher end

**Suspension cells:** These cells need more DNA, please start using the 1.5 x transfection scale and go to 2.0x or 2.5 x.

**Viromer® ONE RED** is fully compatible with our Viromer® RED liquid product, but is easier to handle and significantly reduces variability. Please refer to our cell data base for a complete list of cell types that are amenable to Viromer transfection.



For further information go to [www.viromer-transfection.com/one](http://www.viromer-transfection.com/one) and visit our support pages and cell data bases  
...do not hesitate to contact us!

Version Dec / 16

## Basic transfection protocol

### 96 and 24 well format

Dilute your pDNA/mRNA to 10 ng/ $\mu$ l using sterile H<sub>2</sub>O. Provide a minimum volume of 80  $\mu$ l.

- 1 Pierce the foil and rehydrate each vial with 80  $\mu$ l. Mix swiftly by pipetting up and down and incubate for 15 min at room temperature.
- 2 Transfer complexes to your cells.



transfer volume	5–15 $\mu$ l	25–75 $\mu$ l
pDNA/mRNA per well	50–150 ng	250–750 nl

Incubate cells as usual. Monitor pDNA/mRNA effects 6–48h after transfection. Expression from mRNA can begin as early as 2h.

## Note(s)

For transfections in 12 well or 6 well plates, please rehydrate first and combine units after complexation.

---

---

---

---

---

---

---

---

---

---


---

---

## Transfection Optimization

### 96 well format

Provide 400  $\mu$ l of your pDNA or mRNA (10 ng/ $\mu$ l) in sterile H<sub>2</sub>O. Rehydrate six vials with the volumes indicated below and mix swiftly by pipetting up and down.



Incubate cells as usual. Monitor pDNA/mRNA effects 6–48h after transfection. Expression from mRNA can begin as early as 2h.