

Introduction:

Alzheimer's disease (AD) is a common neurodegenerative disease affecting 35 millions patients worldwide. Microtubule-associated protein tau (Tau) concentration in cerebrospinal fluid (CSF) is a specific and sensitive biomarker for AD. As pathological Tau concentrations are very low (> 0.4 ng/mL), sample preparation must be efficient and Tau quantification must be absolute and accurate to obtain reliable results for preclinical and clinical samples. To achieve this goal, published studies¹ have promoted the use of a full-length stable isotope labelled Tau protein as the gold standard for quantification using LC-MS. Promise Proteomics offers a stable isotope labelled Tau-441 isoform protein standard and this poster presents data on both its characterization and its application in an LC-MS workflow for quantifying Tau in serum samples.


Characterization of the full length Tau-441 protein and its labelled ¹³C, ¹⁵N form

Both labelled and unlabelled full length Tau proteins were produced using an optimized *Escherichia coli* expression system. Produced proteins were analysed using standard quality control procedures: SDS-PAGE; BCA quantitation; intact mass analysis and peptide mapping (SCIEX X500B); isotopic incorporation (SCIEX 6500QTRAP).


Intact mass analysis

Unlabelled and **¹³C,¹⁵N labelled Tau** proteins were separately diluted to a final concentration of 0.1 mg/mL. Analyses were performed on a X500B QToF equipped with a Exion LC System (SCIEX). Intact proteins were separated using a BioZen 3.6 µm Intact C4 150 x 2.1 mm column (Phenomenex).

• Unlabelled Tau at 0.1 mg/mL



• Labelled Tau at 0.1 ng/mL



Observed mass (Da)	Expected mass (Da)	Mass deviation (Da)	Interpretation
45717.22	45849.91	132.69	Unlabelled tau 1-441 without Methionine

Observed mass (Da)	Expected mass (Da)	Mass deviation (Da)	Interpretation
46202.56	46338.35	135.79	Labelled tau 1-441 without Methionine

Peptide Mapping analysis and isotope incorporation

2 µg of unlabelled and labelled Tau protein were digested separately with trypsin. Digested samples were analyzed on a X500B QToF and 6500QTrap, both equipped with a Exion LC System (SCIEX). Peptides were separated using a BioZen™ 2.6 µm Peptide XB-C18, LC Column 100 x 2.1 mm (Phenomenex).

- Percentage coverage for labelled and unlabelled Tau : 62.4 % (275/441)

MAEPRQEFVEMEDHAGTYGLGDRKQGGYTMHQDQEGDGLKESPLQTPTEGSEEPGSETSDAKSTPTAEDVTAP
LVDEGAPGKQAAAQPHTEIPEGTTAEAGIGDTPSLEDAAGHVTAQRMVSKSDGTGSDDKKAKGADGKTIATPRGAA
PPGQKQANATRIPAKTPAPKTPSSGEPKSGDRSGYSSPGSPGTPGSRSTPSPPTREPCKVAVVTRPPKSPSSAKS
RLQAPVPMPLKKNVSKIGSTENLKHQPGGGKQVIINKKLDLSNVQSKCGSKDNKIHVPGGGVQIVYKPVDLKVTSKC
GSLNIHHPGGGQVEVKSEKLDKDRVQSKIGSLDNITHVPGGGKIKIETHKLTFRENAKAKTDHGAIEIVYKSPVSGDTS
PRHLSNVSTGSIDMVDSPQLATLADEVASLAKQGL

Detected peptides were underlined.

- ¹³C, ¹⁵N isotopic incorporation yield for labelled Tau : > 99%

Conclusions:

- Both Tau standards produced were pure.
- N-terminal Methionine was removed on both proteins
- Percentage coverage was high at 62.4%
- Isotopic incorporation for labelled Tau was > 99%.

Quantification of Tau protein

Seven calibrators were prepared between 0.1 - 10 ng/mL in 0.5% serum and spiked with 0.3 ng/mL of labelled Tau 1-441. As shown in Figure 1, Tau protein was extracted from artificial CSF (0.5% serum) by protein precipitation using perchloric acid 70% followed by solid phase extraction of the supernatant with StrataX Polymeric Reverse Phase. After trypsin digestion, Tau peptides were oxidized with H₂O₂ and analysed on Exion LC- QTRAP (Sciex). Tau peptides (Table 1) were separated on BioZen™ 2.6 µm Peptide XB-C18, LC Column 100 x 2.1 mm (Phenomenex). 3 analytical injections were performed.

Figure 1




Table 1


Peptide	Precursor ion (m/z)	Charge	Retention time (min)	Collision energy (V)	Fragments
TPSLTPPTR.light	533.798	2+	8.2	28.1	Y ₈ ⁺ , Y ₇ ⁺ , Y ₆ ⁺
TPSLTPPTR.heavy	538.802	2+	8.2	28.1	Y ₈ ⁺ , Y ₇ ⁺ , Y ₆ ⁺
VQIINK.light	357.729	2+	7.3	21.7	Y ₄ ⁺ , Y ₃ ⁺
VQIINK.heavy	361.736	2+	7.3	21.7	Y ₄ ⁺ , Y ₃ ⁺
IGSLDNITHVPGGGNK.light	526.946	3+	8.6	26.3	Y ₆ ⁺ , Y ₁₅ ²⁺ , Y ₁₄ ²⁺
IGSLDNITHVPGGGNK.heavy	529.617	3+	8.6	26.3	Y ₆ ⁺ , Y ₁₅ ²⁺ , Y ₁₄ ²⁺
SPVSGDTSR.light	551.28	2+	7	28.7	Y ₈ ⁺ , Y ₇ ⁺
SPVSGDTSR.heavy	556.284	2+	7	28.7	Y ₈ ⁺ , Y ₇ ⁺

Results:

Tau peptides TPSLTPPTR, VQIINK, SPVSGDTSR and IGSLDNITHVPGGGNK gave consistent results with each other (Table 2). Calibration curves (Figure 2) were linear from 0.1 to 10 ng/mL and results were repeatable (RSD < 25% for LOQ (0.1 ng/mL) and RSD < 20% for calibrator points).

Conclusion:

Combining Promise Proteomics labelled Tau internal standard with this sample preparation procedure allowed quantification of Tau protein at high sensitivity and this can be applied to pre-clinical or clinical samples.



Chromatogram of IGSLDNITHVPGGGNK (red) and its labelled form (blue) for G2 calibrator (0.5 ng/mL)

Figure 2




Table 2

Tau concentration (ng/mL)	[(Area light/area heavy) * heavy concentration] mean				CV (%)
	pep_TPSLTPPTR	pep_VQIINK	pep_SPV...SPR	pep_IGS...NK	
0.1	0.05	0.06	0.03	0.04	24.60
0.5	0.11	0.10	0.10	0.10	4.87
1	0.19	0.15	0.19	0.19	10.03
2.5	0.46	0.41	0.44	0.38	8.67
5	0.80	0.71	0.73	0.69	6.17
10	1.58	1.54	1.58	1.57	1.17

Authors

Guillaume Picard, Marina Iannello, Pauline Bros, Ravindra Chaudhari, Dorothee Lebert