Whitepaper

“BelDIA PRE”

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Quellen, Publikationen

1. Pfützner A, Standl E, Hohberg C, Konrad T, Strotmann HJ, Lübben G, Langenfeld MR, Schulze J, Forst T. IRIS II study: intact proinsulin is confirmed as a highly specific indicator for insulin resistance in a large cross-sectional study design. *Diabetes Technol Ther.* 2005;7(3):478–486.
2. Pfützner A, Pfützner AH, Larbig M, Forst T. Role of intact proinsulin in diagnosis and treatment of type 2 diabetes mellitus. *Diabetes Technol Ther.* 2004;6(3):405–412.
3. Pfützner A, Pfützner AH, Stridde E, Huppertz E, Reimer T, Derwahl M, Forst T, Petrak F. Insulin resistance and b-cell-dysfunction in insufficiently controlled type 2 diabetes: the SETT2D Trial. *Diabetes Stoffwechsel Herz.* 2007;2:91–97.
4. Pfützner A, Kunt T, Hohberg C, Mondok A, Pahler S, Konrad T, Lübben G, Forst T. Fasting intact proinsulin is a highly specific predictor of insulin resistance in type 2 diabetes. *Diabetes Care.* 2004;27(3):682–687.
5. Pfützner A, Kann PH, Pfützner AH, Kunt T, Larbig M, Weber MM, Forst T. Intact and total proinsulin - new aspects for diagnosis and treatment of type 2 diabetes mellitus and insulin resistance. *Clin Lab.* 2004;50(9-10):567–573.
6. Galloway JA, Hooper SA, Spradlin CT, Howey DC, Frank BH, Bowsher RR, Anderson JH. Biosynthetic human proinsulin. Review of chemistry, in vitro and in vivo receptor binding, animal and human pharmacology studies, and clinical trial experience. *Diabetes Care.* 1992;15(5):666–692.
7. Pfützner A, Kunt T, Langenfeld M, Löbig M, Knesovic M, Forst T. Clinical and laboratory evaluation of specific chemiluminescence assays for intact and total proinsulin. *Clin Chem Lab Med.* 2003;41(9):1234–1238.
8. Pfützner A, Pfützner AH, Kann PH, Stute R, Löbig M, Yang JW, Mistry J, Forst T. Clinical and laboratory evaluation of a new specific ELISA for intact proinsulin. *Clin Lab.* 2005;51(5-6):243–249.
9. Siebenhaar R, Weise A, Safinowski M, Reisinger K, Musholt PB, Reimer T, Pfützner A, Forst T. Clinical and laboratory evaluation of a new specific ELISA for intact proinsulin. *Diabetes Stoffwechsel Herz.* 2008; 2:275–281.
10. Pfützner A., Pfützner AH., Kann PH, Burgard G. Clinical and laboratory evaluation of a new specific point-of-care test for intact proinsulin. *J. Diabetes Sci. Technol.* 2017;11(2):278-283
11. Pfützner A, Hermanns I, Ramljak S, Demircik F, Pfützner AH, Kann PH, Weber MM. Elevated intact proinsulin during an oral glucose challenge indicate progressive β -cell dysfunction and may be predictive for development of type 2 diabetes. *J. Diabetes Sci. Technol.* 2015; 9(6):1307-1312
12. Nijpels G, Popp-Snijders C, Kostense PJ, Bouter LM, Heine RJ. Fasting proinsulin and 2-h post-load glucose levels predict the conversion to NIDDM in subjects with impaired glucose tolerance: the Hoorn Study. *Diabetologia.* 1996;39(1):113–118.
13. Alsema M, Dekker JM, Nijpels G, Stehouwer CD, Bouter LM, Heine RJ; Hoorn Study. Proinsulin concentration is an independent predictor of all-cause and cardiovascular mortality: an 11-year follow-up of the Hoorn Study. *Diabetes Care.* 2005;28(4):860–865.
14. Nordt TK, Bode C, Sobel BE. Stimulation in vivo of expression of intra-abdominal adipose tissue plasminogen activator inhibitor type I by proinsulin. *Diabetologia.* 2001;44:1121-1124
15. Zethelius B, Byberg L, Hales CN, Lithell H, Berne C. Proinsulin is an independent predictor of coronary heart disease: report from a 27-year follow-up study. *Circulation.* 2002;105:2153-2158
16. Bavenholm P, Proudler A, Tornvall P, Godsland I, Landou C, de Faire U, Hamsten A. Insulin, intact and split proinsulin, and coronary artery disease in young men. *Circulation.* 1995;92(6):1422–1429
17. Wareham NJ, Byrne CD, Hales CN. Role of insulin and proinsulin in diabetic vascular disease. *Metabolism.* 1995;44(10 Suppl 4):76–82.
18. Lindahl B, Dinesen B, Eliasson M, Røder M, Jansson JH, Huhtasaari F, Hallmans G. High proinsulin concentration precedes acute myocardial infarction in a nondiabetic population. *Metabolism.* 1999;48(9):1197–1202.
19. Lindahl B, Dinesen B, Eliasson M, Røder M, Hallmans G, Stegmayr B. High proinsulin levels precede first-ever stroke in a nondiabetic population. *Stroke.* 2000;31(12):2936–2941.
20. Yudkin JS, May M, Elwood P, Yarnell JW, Greenwood R, Davey Smith G; Caerphilly Study. Concentrations of proinsulin like molecules predict coronary heart disease risk independently of insulin - prospective data from the Caerphilly Study. *Diabetologia.* 2002;45(3):327–336.
21. Katz RJ, Ratner RE, Cohen RM, Eisenhower E, Verme D. Are insulin and proinsulin independent risk markers for premature coronary artery disease *Diabetes.* 1996;45(6):736–741.

22. Røder ME, Porte D, Schwartz RS, Kahn SE. Disproportionately elevated proinsulin levels reflect the degree of impaired β -cell secretory capacity in patients with non-insulin dependent diabetes mellitus. *J Clin Endocrinol Metab.* 1998;83:604-608.
23. Haffner SM, Mykkänen L, Stern MP, Valdez RA, Heisserman JA, Bowsher RR. Relationship of proinsulin and insulin to cardiovascular risk factors in nondiabetic subjects. *Diabetes.* 1993;42:1297-1302.
24. Langenfeld MR, Forst T, Standl E, Strotmann HJ, Lübber G, Pahler S, Kann P, Pfützner A; IRIS II study. IRIS II Study: sensitivity and specificity of intact proinsulin, adiponectin, and the proinsulin/adiponectin ratio as markers for insulin resistance. *Diabetes Technol Ther.* 2004;6(6):836–843. [[PubMed](#)]
25. Pfützner A, Kunt T, Mondok A, Pahler S, Konrad T, Lübber G, Forst T. Fasting intact proinsulin is a highly specific predictor of insulin resistance in type 2 diabetes. *Diabetes Care.* 2004 ;27:682-687.

Contact

PharmACT AG for BelDIA PRE

Prof. Dr. Dr. med. Andreas Pfützner, PHSI Institute, Mainz, Germany.

Mailto: customer@pharmact.eu

//ENDE//