

Medela Healthcare Accessories



MEDICAL VACUUM TECHNOLOGY FOR HEALTHCARE PROFESSIONALS

Proven quality

The Medela Healthcare range of products fulfils international industrial standards: ISO 9001 and ISO 13485 as well as the American Good Manufacturing Practices guidelines (GMP). Our products carry the CE marking in accordance with the directive 93/42/EEC for medical products.

For us, the assurance of a long guarantee is taken for granted.

Symbols



This symbol indicates a single use device. Do not reuse the device.



Autoclavable – for example at 134°/10 minutes



Not autoclavable – cleaning in a cleaning and disinfecting machine recommended

Clario

4



The Clario system is ideal for airway suctioning at home, on the road or in the clinic. Easy to use and highly mobile.

Vacuum settings:

min. – 18 kPa – 135 mmHg
med. – 36 kPa – 270 mmHg
max. – 80 kPa – 600 mmHg

Flow: 15l/min

General Accessories



Ref. 074.0001

Canister set

For collecting secretions when using the Clario pump
– Set for Clario with canister 0.5l, lid with plugs and floater



Ref. 074.0002

SafetyChamber

For vacuum regulation
– SafetyChamber for Clario with vacuum regulator knob and 2 valve discs



Ref. 074.0003

Disposable suction tubing with fingertip, sterile

Attaches suction catheter to pump for suctioning purposes
– PVC
– Length 150 cm, CH 25, sterile



Ref. 074.0006

Vacuum gauge

For measuring vacuum levels of the Clario pump
– Can be plugged into lid of Clario canister
– In kPa and mmHg



Ref. 074.0007

SoftCase carrying bag

For safe transportation of a home care pump and accessories
– With 2 carrying straps for attachment to a wheelchair or as a backpack and handle



Ref. 074.0017

Car connection cable

For operating the Clario AC/DC pump in the car
– 12V DC cable

Clario Toni

Just the right paediatric care – easy to use, effective and highly mobile.



Vacuum settings:

- min. – 60 mmHg – 8 kPa
- med. – 100 mmHg – 13 kPa
- max. – 600 mmHg – 80 kPa

Flow: 15l/min

General Accessories



Ref. 074.0001

Canister set

For collecting secretions when using the Clario pump
 – Set for Clario with canister 0.5l, lid with plugs and floater



Ref. 014.0030

SafetyChamber Clario Toni

For vacuum regulation
 – SafetyChamber for Clario Toni with vacuum regulator knob and 2 valve discs



Ref. 074.0003

Disposable suction tubing with fingertip, sterile

Attaches suction catheter to pump for suctioning purposes
 – PVC
 – Length 150 cm, CH 25, sterile



Ref. 074.0006

Vacuum gauge

For measuring vacuum levels of the Clario pump
 – Can be plugged into lid of Clario canister
 – In kPa and mmHg



Ref. 074.0007

SoftCase carrying bag

For safe transportation of a home care pump and accessories
 – With 2 carrying straps for attachment to a wheelchair or as a backpack and handle



Ref. 074.0017

Car connection cable

For operating the Clario AC/DC pump in the car
 – 12V DC cable

Vario 8, Vario 8 c/i, Vario 18 and Vario 18 c/i

6



The Vario 8 and Vario 8 c/i low-vacuum suction pumps are all-rounders for hospitals and outpatient clinics. They offer versatility, reliability and mobility and are suitable for a wide range of low-vacuum applications. Vario 18 is an efficient high-vacuum medical suction pump. As an all-rounder for hospitals, clinics and home care, it offers reliability, low noise and mobility. The Vario 18 c/i is a medium-vacuum medical suction pump for hospitals and clinics, providing constant and intermittent suction.

	Vacuum settings:	Flow:
Vario 8, Vario 8 c/i:	-9 kPa -68 mmHg	8 l/min
Vario 18 c/i:	-55 kPa -413 mmHg	18 l/min
Vario 18:	-75 kPa -563 mmHg	18 l/min

General Accessories

Can be used together with Medela Fluid Collection Systems (Page 13/14), Tubings (Page 17), Filters (Page 18) and Connectors (Page 19)



Ref. 077.0511

Trolley small

- For transport and use of various accessories
- Nickel-plated with 2 standard rails 10 x 25 mm
 - Mobile with 4 castors (2 of which have brakes)
 - Dimensions: 45 x 41 x 86 cm



Ref. 077.0521

Clamp holder

- For attaching suction jars to standard rails
- Polyamid



Ref. 077.0821

Vario rail holder 10 x 25 mm

- For attaching Vario to standard rails
- Aluminum holder



Ref. 077.0823

Vario universal holder

- For attaching the Vario to rails
- Holder for Vario can be rotated to 90°



Ref. 077.0831

Carrying bag

- For safe transportation of the Vario and accessories
- Sturdy carrying case with shoulder strap and pockets



Ref. 077.0832

Car connection cable

- For operating the Vario AC/DC pump in the car
- 12V DC cable

Basic



The Basic provides reliable and simple suction for multiple applications in medical practices, hospital wards, operating theatres and emergency departments.

Vacuum settings:

max. -90 kPa - 675 mmHg

Flow: 30l/min

General Accessories

Can be used together with Medela Vacuum Assisted Delivery System (Page 15/16), Fluid Collection Systems (Page 13/14), Tubings (Page 17), Filters (Page 18) and Connectors (Page 19)



Ref. 017.0015

Apgar Timer

- For reliable notification to perform Apgar tests
- Battery operated
 - Acoustic signal after 1, 5 and 10 minutes
 - Can be attached to standard rails



Ref. 077.0521

Clamp holder

- For attaching suction jars to standard rails
- Polyamid



Ref. 077.0532 (480 mm)

Ref. 077.0531 (280 mm)

Quiver with holder

- For storing catheters or suctioning instruments during an operation
- 2 parts. Holder (PC) and quiver (PP) can be attached directly to standard rails
 - Length: 280 mm or 480 mm



Ref. 077.0723

Foot on/off switch


- For turning the Basic on/off
- Cable length 3.5 m



Ref. 077.0731

Foot vacuum regulator

- For regulating the vacuum of the Basic
- Can be attached to standard Medela suction jars
 - With hook and clamp holder for easy attachment to standard rails when not in use



Ref. 077.0751 (l)

Ref. 077.0752 (r)

Tubing holder

- For attaching tubing to standard rails
- With secure attachment screw
 - For either right or left-sided attachment

Dominant Flex

8



Designed for use in hospitals, clinics and medical practices, the Dominant Flex is the optimal medical suction pump for all suctioning needs.

Vacuum settings:

max. -95 kPa - 713 mmHg

Flow: 40 l/min
50 l/min
60 l/min

General Accessories

Can be used together with Medela Vacuum Assisted Delivery System (Page 15/16), Fluid Collection Systems (Page 13/14), Tubings (Page 17), Filters (Page 18) and Connectors (Page 19)



Ref. 077.0731

Foot vacuum regulator

For regulating the vacuum of the Dominant Flex

- Can be attached to standard Medela suction jars
- With hook and clamp holder for easy attachment to standard rails when not in use



Ref. 077.0723

Foot on/off switch

For turning the Dominant Flex on/off

- Cable length 3.5 m



Ref. 077.0521

Clamp holder

For attaching suction jars to standard rails

- Polyamid




Ref. 077.0532 (480 mm)

Ref. 077.0531 (280 mm)

Quiver with holder

For storing catheters or suctioning instruments during an operation

- 2 parts. Holder (PC) and quiver (PP) can be attached directly to standard rails
- Length: 280 mm or 480 mm



Ref. 077.0751 (l)

Ref. 077.0752 (r)

Tubing holder

For attaching tubing to standard rails

- With secure attachment screw
- For either right or left-sided attachment



Ref. 017.0015

Apgar Timer

For reliable notification to perform Apgar tests

- Battery operated
- Acoustic signal after 1, 5 and 10 minutes
- Can be attached to standard rails

Thopaz+™



Medela's Thopaz+ takes chest drainage therapy to a new level of care. Unlike analogue systems, it reliably regulates the applied pressure at the patients chest and digitally monitors crucial therapy indicators.

Vacuum settings:

max. - 10 kPa (+/- 15%)

Flow: max. 5l/min (+/- 10%)

Thopaz+ Disposables

Together with a canister and a tubing Thopaz+ is a complete cardi thoracic drainage system.



Ref. 079.0011 (0.3l)
 Ref. 079.0016 (0.8l)
 Ref. 079.0018 (2.0l)

Canister without solidifier, sterile

Standard canisters for the collection of fluids



Ref. 079.0012 (0.3l)
 Ref. 079.0017 (0.8l)
 Ref. 079.0019 (2.0l)

Canister with solidifier, sterile

Provides a hygienic and a cost effective disposal



Ref. 079.0021 (single)
 Ref. 079.0022 (double)

Tubings, sterile

With different connectors to the drain
 Length 1.5 m



Ref. 079.0024 (double)
 Ref. 079.0023 (single)

Tubings with small connectors, sterile

For the connection with small drains
 Length 1.5 m



Ref. 079.0025 (single)
 Ref. 079.0026 (double)

Tubings with large connectors, sterile

For the connection with large drains
 Length 1.5 m



Ref. 079.0039

Sealing cap, sterile

For closing an unused connection

Thopaz+ Accessories

The optimally designed accessories can be used for charging the battery, carrying and attaching Thopaz+ on rails, beds and IV poles.



Ref. 087.0059

Wall mains adapter

For charging the battery

With different plugs (Euro-UK-US-AUS/NZ)



Ref. 079.0037

Docking station

For holding Thopaz+ and charging the battery

Can be attached to standard rails, beds* and IV poles*

*using a holder with standard rail



Ref. 079.0038

Adapter Docking station for Thopaz+ canister 2 liter

To secure the 2l canister in position



Ref. 079.0036

Holder with standard rail

For the easy attachment of Thopaz on rails, beds, IV poles

Can be fixed in a vertical or horizontal position



Ref. 079.0040

Universal holder with flexible standard rail

For the easy attachment of Thopaz+ on rails, beds, IV poles

Stepless adjustment of the inclination of the rail (360°)



Ref. 079.0034

USB Cable

For the connection of Thopaz+ with the PC

Type: 2.0 A/M- Mini B/M



Ref. 079.1006

Carrying strap

For carrying Thopaz+

Invia Liberty, Invia Motion and Invia Dressing Kits*



Medela offers efficient Negative Pressure Wound Therapy solutions, simplifying the delivery and management of NPWT – making it easier for healthcare professionals to help their patients heal.

Medela offers a full range of accessories for its NPWT solutions.

Invia Liberty Accessories



Ref. 087.0012 (300ml)
Ref. 087.0017 (800ml)

Canisters

- Two available sizes (300ml and 800ml) for maximum flexibility.
- Built-in bacteria filter.
- Self-sealing filter protects pump from liquids.
- Built-in solidifier for convenient disposal.
- Sterile, for single use.



Ref. 087.0025

Tubing with Quick-connector

- Double lumen tubing for active exudate management and pressure control.
- Tubing with Quick-connector for easy and secure attachment



Ref. 087.0031

Carrying case

- Designed for discreet and convenient pump transportation.
- Easy opening for canister changes.
- Accessible display even when the bag is closed.



Ref. 087.0059

Power Supply (EU)

- Replacement available separately.
- Only use Medela provided power supply with Invia Liberty devices.



Ref. 079.0037

Docking Station

- Convenient placement and removal for easy pump charging.
- Can be placed on level surface or clamped to a rail.

* Medela NPWT solutions are currently not available in all markets

Invia Motion accessories



Ref. 087.4016

Canister/Tubing Set 0.15 l

- Disposable canister/tubing Set.
- Available size: 150ml.
- Sterile, for single use.



Ref. 087.4003

Carrying Case

- The Invia Motion System is intended to be worn in the carrying case during operation.
- The carrying case is provided with the Invia Motion System.
- Shoulder/belt strap included.




Ref. 200.8604

Power Supply

- Power supply provided in the Invia Motion starter kit; replacement available separately.
- Only use Medela provided power supply with Invia Motion devices.

Invia Dressing Kits



Ref. 087.0023


Invia Y connector



Ref. 087.0028

Invia FitPad


Single sterile Fitpad suction interface



Ref. 087.0030

Invia Transparent Film

Single sterile Transparent films (26cm x 32cm)



- Ref. 087.6221 (small, 3 pcs)
- Ref. 087.6225 (small, 15 pcs)
- Ref. 087.6222 (medium, 3 pcs)
- Ref. 087.6226 (medium, 15 pcs)
- Ref. 087.6223 (large, 3 pcs)
- Ref. 087.6227 (large, 15 pcs)
- Ref. 087.6224 (X-large, 3 pcs)

Invia Foam Dressing Kit with FitPad



- Ref. 087.7170 (medium, 3 pcs)
- Ref. 087.7172 (medium, 15 pcs)
- Ref. 087.7171 (large, 3 pcs)

Invia Gauze Dressing Kit with FitPad

Medela DCS



Medela's Disposable Collection System supports the safe collection and disposal of suctioned fluids.

Disposable Collection System



Ref. 077.0082 (1.5l)
Ref. 077.0085 (2.5l)

Suction jar PC for disposable liners

Reusable jars for holding disposable liners
 – Made of unbreakable and easy to clean polycarbonate
 – Sizes: 1.5 and 2.5 liter



Ref. 077.0083 (1.5l)
Ref. 077.0086 (2.5l)

Disposable suction liners

Disposable solution for collecting secretions
 – 2 layered foil with integrated overflow protection and lid
 – Sizes: 1.5 and 2.5 liter



Ref. 077.0084 (1.5l)
Ref. 077.0087 (2.5l)

Disposable suction liners with solidifier

Disposable solution for collecting secretions
 – 2 layered foil with integrated overflow protection and lid
 – Contains solidifier for leak-proof transport and disposal
 – Sizes: 1.5 and 2.5 liter



Ref. 077.0094

Disposable specimen cup

For collecting tissue samples
 – Can be attached to patient connection of lid
 – Made of polypropylene
 – Package of 10 pcs



Ref. 077.0188 A
Ref. 077.0190 B

Trolley

For use with Medela Fluid Collection Systems
 – Holds up to 5 suction jars/liners
 – Mobile with 5 castors (with brakes)
 – Dimensions: 85 x 60 x 60 cm

Medela Reusable Collection System



The Medela reusable collection system, made out of durable, high quality components, facilitates the collection of suctioned fluids.

Reusable Collection System



Ref. 077.0125 (0.25 l)
 Ref. 077.0155 (0.5 l)
 Ref. 077.0110 (1 l)
 Ref. 077.0120 (2 l)
 Ref. 077.0130 (3 l)
 Ref. 077.0150 (5 l)

Suction jar, PSU

Reusable jar for collecting suctioned secretions

- Made of unbreakable polysulfone with easy-to-read white graduation
- Sizes 0.25, 0.5, 1, 2, 3 and 5 liter



Ref. 077.0420 (Ø 6-10 mm)
 Ref. 077.0430 (Ø 10-14 mm)

Lid large

Lid for reusable jars

- Lid with seal for use with the suction jars 1, 2, 3 and 5 liter with overflow protection device
- Either with Ø 6 – 10mm or 10 – 14mm conical patient connection



Ref. 077.0440

Lid small

Lid for reusable jars (0.25 and 0.5 l)

- Lid with seal for use with the suction jars 0.25 and 0.5 liter
- Ø 6 – 10mm conical patient connection and overflow protection device



Ref. 077.0702 (2 l)
 Ref. 077.0701 (1 l)
 Ref. 077.0703 (3 l)
 Ref. 077.0705 (5 l)

Reusable set, suction jar PSU

Complete set reusable for collecting secretions consists of:

- Suction jar, PSU, lid with conical patient connection Ø 6 – 10mm and overflow protection device, plastic clamp holder, disposable tubing 180 cm with fingertip
- Sizes: 1, 2, 3 and 5 liter



Ref. 077.0562 (Ø 6-10 mm)
 Ref. 077.0564 (Ø 10-14 mm)

Disposable specimen cup

For collecting tissue samples


- Can be attached to patient connection of lid
- Made of polypropylene
- Package of 10 pcs

Medela VAD



The Medela Vacuum Assisted Delivery System is an individual and efficient system for safe interventions. Together with the Medela Basic suction pump and the choice of rigid or soft, reusable or disposable suction cups and a wide range of accessories, it delivers safe, individualised assistance.

Reusable Suction Cups



Ref. 077.0177 (Ø 40 mm)
 Ref. 077.0176 (Ø 50 mm)
 Ref. 077.0175 (Ø 60 mm)

Malmström Cup

- Autoclavable stainless steel
- Used when the fetus is in the occiputo-anterior position
- Cup design allows for good seal to help minimise pop-offs
- Sizes: Ø 40, 50 and 60 mm



Ref. 077.0173 (Ø 40 mm)
 Ref. 077.0172 (Ø 50 mm)
 Ref. 077.0171 (Ø 60 mm)

Bird Cup

- Autoclavable stainless steel
- Based on the design of the Malmström Cup, with separate suction and traction ports
- Sizes: Ø 40, 50 or 60 mm



Ref. 077.0174

Bird Posterior Cup

- Autoclavable stainless steel
- Especially designed for use when the fetus is in the occiputo-posterior position
- Vacuum port on side of suction cup and low profile to ensure easier insertion
- Size: Ø 50 mm



Ref. 077.0761 (Ø 50 mm) blue
 Ref. 077.0762 (Ø 60 mm) yellow

Siic Cup

- Soft suction cup for vacuum assisted delivery
- With silicone tubing, 1.5 m and 1 coupling piece
- Sizes: Ø 50 and 60 mm



Ref. 077.0078

Siic Cup with trumpet valve

- Attach gently to fetal scalp
- One-piece design for easy handling and use
- Smooth surface with longitudinal ridge makes it possible to observe rotation
- Size: Ø 60 mm




Ref. 077.0080

CaesAid Cup for Caesarean sections

- Autoclavable medical grade silicone
- Gentle support for delivering fetal head during Caesarean sections
- Can be attached to the pump used for suctioning during the operation.
- Size: Ø 60 mm


Disposable Suction Cups



Ref. 077.0791

Single Use Silc Cup

- Attaches gently to fetal scalp
- Smooth surface with longitudinal ridge makes it possible to observe rotation
- Sterile with 150 cm tubing
- Size: Ø 60 mm



Ref. 077.0792


Single Use Bird Cup

- Polycarbonate version of the reusable Bird Cup
- Traction handle can be tilted for use when the fetus is in the occiputo-anterior or posterior position
- Sterile with 150 cm tubing
- Size: Ø 50 mm

Tubings


For hygienic transportation of secretions, reliable tubings complete the high quality system. Depending on the application, Medela offers various materials, diameters and connection pieces with the tubing.

Silicone tubings



Ref. 077.0901 (1 m)
 Ref. 077.0902 (2 m)
 Ref. 077.0905 (5 m)
 Ref. 077.0900 (25 m)

Silicone tubings
 – Ø 7 x 12 mm
 – Lengths 1, 2, 5, or 25 m



Ref. 077.0912 (60 cm)
 Ref. 077.0913 (100 cm)

Silicone tubings with 1 coupling piece
 – Ø 7 x 12 mm
 – Lengths 60 or 100 cm



Ref. 077.0922 (60 cm)
 Ref. 077.0921 (30 cm)

Silicone tubings with 2 coupling pieces
 – Ø 7 x 12 mm
 – Lengths 30 or 60 cm



Ref. 077.0931
 Ref. 077.0185 (60 cm)

Silicone tubing with 1 coupling piece and 1 coupling piece 90° angle
 – Ø 7 x 12 mm
 – Length 25 cm



Ref. 077.0970

Silicone tubing
 – Ø 5 x 10 mm
 – Length 25 m

Disposable tubing



Ref. 077.0951

Disposable tubing with fingertip
 – PVC, sterile
 – Ø 6 x 9 mm
 – Length 180 cm



Ref. 077.0965
 Ref. 077.0952
 Ref. 077.0191

Disposable tubing with funnel connectors, Ø 5.5 x 8.5 mm
 – with two funnel connections, 150 cm
 – with one funnel connection (non-sterile), 180 cm
 – with one funnel connection, 300 cm



Ref. 077.0184

Disposable tubing with universal connectors
 – for use with Single Use Bird Cup or reusable Bird Cups
 – PVC, sterile
 – Ø 5.5 x 8.3 mm
 – Length 150 cm

Filters

18

Customise your suction system according to your hygienic needs with filters. They are used to prevent bacteria growth, the development of odours or for overflow protection.

Disposable Filters



Ref. 077.0092

Disposable overflow protection/bacteria/odour filter

Neutralizes odour, prevents an overflow and protects from bacteria

- With Medela connections
- Package of 10 pcs



Ref. 077.0571 (Medela connections)
Ref. 077.0572 (conical connections)
Ref. 077.0577 (Medela and conical connections)

Disposable overflow protection / bacteria filter

Prevents an overflow and protects from bacteria

- Package of 10 pcs
- Choice of Medela or conical connections



Ref. 077.0573 (Medela connections)
Ref. 077.0574 (conical connections)

Disposable bacteria filter

Protects pump and environment from bacteria

- Package of 10 pcs
- Choice of Medela or conical connections



Ref. 077.0575 (Medela connections)
Ref. 077.0576 (conical connections)

Disposable bacteria/odour filter

Filters and neutralizes odour and protects from bacteria

- Package of 10 pcs
- Choice of Medela or conical connections

Connections

Connecting tubings to each other or to accessories often requires connection pieces. These have been carefully developed to fulfill the requirements for a tight connection and an easy fit.

Reusable and disposable connections



Ref. 077.1015 (Ø 6–10 mm) **Angle piece**
 Ref. 077.1016 (Ø 10–14 mm) Connects lid/suction liners to patient tubing

- and prevents bending
- Polypropylene
- Package of 10 pcs
- Sizes: Ø 6–10 or 10–14 mm



Ref. 077.1017 **Double conical coupling piece**

- Connects tubings to each other
- Polypropylene
- Package of 10 pcs



Ref. 077.1018 **Coupling piece**

- Connects tubing to Medela connections (lid, change-over-valve, etc.)
- Polyoxymethylen with 2 blue silicone o-rings
- Package of 10 pcs



Ref. 077.1019 **Coupling piece 90°**

- Connects tubing to Medela connections (lid, change-over-valve, etc.)
- Polyoxymethylen with 2 blue silicone o-rings
- Package of 10 pcs



Ref. 077.1020 **Y-piece**
 For attaching 2 tubings to the same vacuum source

- Polypropylene
- Package of 10 pcs

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.



Medical Vacuum Technology for Healthcare Professionals

**Please contact us or your local Medela
representative for details.**



Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact:

Australia

Medela Pty Ltd
Medical Technology
3 Arco Lane
Heatherton, Vic 3202
Australia
Phone +61 (0) 3 9552 8600
Fax +61 (0) 3 9552 8699
contact@medela.com.au
www.medela.com.au

Canada

Medela Inc.
4160 Sladeview Cres., Unit #8
Mississauga, Ontario
L5L 0A1
Canada
Toll Free 800 435 8316
Toll Free Fax 800 995 7867
info@medela.ca
www.medela.ca

India

Medela India private limited
c/o Vatika Business Park
First floor, tower 2,
Sohna Road, Sec-49,
Gurgaon 122 002
India
Phone +91 124 4416999
Fax +91 124 4416990
info@medela.in
www.medela.in

UK

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
UK
Phone +44 (0)161 776 0400
Fax +44 (0)161 776 0444
info@medela.co.uk
www.medela.co.uk

USA

Medela Inc.
1101 Corporate Drive
McHenry, IL 60050
USA
Phone +1 877 735 1626
Fax +1 815 307 8942
info-healthcare@medela.com
www.medela-healthcare.us