

Rodneys aligner seater TM

- The Rodney's Aligner Seater with rod helps aligner trays fit closer to the teeth to ensure it reaches its preset position and orthodontic force.
- The unique plastic rod is designed for easier control and holds the soft cylinder in place.
- The rod is more convenient and hygienic for the patient who no longer need to put their fingers in the mouth.
- The rod can be easily removed.
- Using Rodney's regularly will ensure a proper fit for the aligners. We recommend using Rodney's 3 to 5 times per day for about 2 to 3 minutes and each time a new aliner is used.
The goal is to ensure the aligners fit correctly and tightly with no gaps in between.
- Rodney's are helpful for teeth with less surface area, for example, teeth with smaller crowns.
- They are also helpful in locking in attachments to the aligners.
- Made of food class silicone.
- Comes in 4 different flavors, mint, lemon, watermelon and unscented.



Rodneys aligner seater	
Description	Part #
Mint scented	RAS-001
Lemon scented	RAS-002
Watermelon scented	RAS-003
Unscented	RAS-004

