

VITRUMAN™

Ultimate Aspirator B89025



b BLACK & BLACK
S U R G I C A L

VITRUVIAN™ Ultimate Aspirator B89025

The Black & Black Vitruvian™ Ultimate aspirator is one of the most powerful, quiet and versatile liposuction aspirators on the market. Powered by two oil-less, rocking piston vacuum pumps, this machine provides excellent suction power along with quiet operation. Its dual voltage motor can be used with 120V/230V electrical service.

Our Ultimate aspirator is carefully designed for optimal performance, with an upright stance and shopping-cart style handle for ease of mobility. Bright LED vacuum



gauges are easy to read. The collection canisters are mounted for ease of viewing during surgery. The single knob control can be adjusted 0 to 29.5 in of Hg vacuum.

While the Vitruvian Ultimate aspirator can be used with all brands of collection canisters, for best performance use Black & Black liposuction canisters (B89140).

- Large 3" hospital grade wheels
- Dual pneumatic foot-switches
- Dual voltage motors 120V/230V
- Bright, easy to read vacuum gauges
- Dual motors for optimal performance
- Quiet operation (50 dB)
- Easy to see collection canisters
- Extendable pole to hold fluid bags
- Convenient storage basket

B89025	VITRUVIAN™ ULTIMATE ASPIRATOR, DUAL MOTOR
B89140	VITRUVIAN™ LIPOSUCTION CANISTERS; SINGLE USE, NON-STERILE, 2000CC, 12/CASE
B89104	ASPIRATOR FILTER, .3 MICRON FILTER, HEPA FOR ASPIRATOR B 89025
B89200	HIGH FLOW ASPIRATION TUBING, 10' - BOX OF 10

#25-3014 Rev:B 2017.02



*Infiltrator pump not included



5238 Royal Woods Parkway, Suite 170
 Tucker, GA 30084
 770.414.4880 phone
 1.877.252.2517 toll-free
 770.414.4879 fax
 info@blackandblacksurgical.com
 www.blackandblacksurgical.com