

# Single-spot LED Matrix Light

with electronic dimmer



Quality of Inspection

LED illuminator containing 18 ultra-bright white LED lamps mounted on gooseneck arm for various oblique incident light applications. Equipped with a universal clamp compatible with all Inspectis microscope Stand models.



Optimised for illumination of large and shiny surfaces like electronics PCBs and metal parts.

Spot-light is powered directly via AUX output of Inspectis Digital Microscopes for excellent cable management. ON-OFF switch and brightness control is integrated. Extra power supply for compatibility to other devices are provided as option.

## Specification:

Number of light elements	18 pcs Ø5 mm LED lamps
Colour temperature	8600 K
Size of LED panel	78x32.5x23 mm
Power	12VDC, 300mA via AUX or External power supply
Length of gooseneck arm	120mm
Total weight	280 g. max

## Ordering Number/Name

HD-125      Inspectis Single-Spot LED Matrix Light  
HD-024-S      External Power Supply (optional)

Reseller:

