

# Can a dialyzer expedite extracorporeal gas exchange?

- **F. Mouzakis<sup>1</sup>, A. Kashefi<sup>1</sup>, S. Rütten<sup>2</sup>, J. Spillner<sup>3</sup>, K. Mottaghy<sup>1</sup>, F. Hima<sup>3</sup>**
- <sup>1</sup> ECC Lab Institute for Physiology, University Hospital Medical Faculty, RWTH Aachen University, Aachen, Germany
- <sup>2</sup> Electron Microscopy Facility, Institute of Pathology, RWTH University Aachen, Germany
- <sup>3</sup> Department of Thoracic Surgery, University Hospital Medical Faculty, RWTH Aachen University, Aachen, Germany

# Oxygenator: a luxury single use product



# Status Quo

- Standard hollow fibers for gas exchange
  - Polypropylene (PP): CPB oriented
  - Polymethylpentene (PMP): gold standard for ECMO
    - Fiber mats often develop **microcracks** near the thread knots
- Alternative materials
  - Silicone flat sheets (e.g. Kolobow): still occasionally used in pediatrics
  - Whole silicone fibers: ~35-55  $\mu\text{m}$  wall thickness – reduced gas diffusion rate

# Silicone-coated Fibers

- Microporous hollow fiber membrane as scaffold

> [Ann Thorac Surg.](#) 1997 Jun;63(6):1730-6. doi: 10.1016/s0003-4975(97)00119-7.

## **Silicone-coated polypropylene hollow-fiber oxygenator: experimental evaluation and preliminary clinical use**

T Shimono <sup>1</sup>, Y Shomura, I Hioki, A Shimamoto, H Tenpaku, Y Maze, K Onoda, M Takao, H Shimpo, I Yada



Affiliations + expand

PMID: 9205175 DOI: 10.1016/s0003-4975(97)00119-7

# Silicone-coated Fibers

- Highly porous scaffold materials
- Defined silicone coating thickness
  - Penetration depth vs. Silicone film
- Intra-/extracapillary coating
  - Affordable dialysis membranes (PSU/PES)
  - Conventional dialyzers
    - Suitable for plug & play applications, e.g. CO<sub>2</sub> elimination (ECCO<sub>2</sub>R) - Hypercapnia

# The manufacturing process




Coating with red-dyed silicone for enhanced visibility

# Validation techniques


- Leakage control
  - Gas pressure build-up
- Pressure drop monitoring
  - Gas
  - Liquid
- Scanning electron microscopy (SEM)

# SEM Images






# SEM Images




# Gas exchange performance

- Silicone Coated Dialyzer
  - $A_{\text{eff}}$ : 0.7 / 1.4 / 1.8 m<sup>2</sup>
- PMP oxygenator
  - $A_{\text{eff}}$ : 1.35 m<sup>2</sup>
- Test conditions
  - ISO 7199:2016 - Gas:Blood 1:1
  - COPD: pCO<sub>2</sub> ≥ 80 mmHg, sO<sub>2</sub> ≤ 50%, Gas:Blood 1:1 - 15:1




# O<sub>2</sub> Transfer Rate

ISO7199 Gas:Blood 1:1




# CO<sub>2</sub> Transfer Rate

ISO7199 Gas:Blood 1:1



# Specific CTR

Hypercapnia/COPD  $Q_B$ : 0.35 l/min



# Epilogue

- Silicone-coated fibers/dialyzers
  - efficient & durable coating
  - bio-/hemocompatible surface
  - high gas permeability (comparable-superior to PP/PMP)
  - cost-effective & expeditious production
  - operated by dialysis machines for Low-Flow applications (e.g. ECCO<sub>2</sub>R)
- Prospective implementation
  - Applicability in High-Flow therapies (CPB, ECMO etc.)
  - In vivo animal experiments

рахмат  
 barpnanaa  
 спасибо  
 faafetai lava  
 mersi  
 kinua barka  
 welalin tack  
 ngiyabonga  
**danke**  
 謝謝  
 spas  
 tesekkür ederim  
 mahalo  
 tapadh leat  
 vinaka  
 спасиби  
 biagodaram  
 dank je  
 misaotra  
 matondo  
 paldies grazzi  
 nanni  
 kiitos  
 dankie  
 dhanyavad  
 maururu  
 koshonom  
**thank you**  
 akun  
 dankon  
 aciu  
 go raibh maith agat  
 mamnun  
 mochchakkeram  
 x b a n a  
 asante  
 manana  
 obligada  
 tenki  
 chokrane  
 murakoze  
 enkosi  
 bayarlalaa  
 grazie  
 hvata  
 chnorakaloutioun  
 gratias age  
 gracias  
 suipay  
**obrigado**  
 sobodi  
 dekuji  
 sagulun  
 sukriya  
 kop khun krap  
 taiku  
 go raibh maith agat  
 dakujem  
 trugarez  
 arigato  
 takk  
 dhanyavadagalu  
 shukriya  
 merce  
 мерси  
 merci  
 xoxoxo  
 diolch  
 tanemirt  
 rahmet  
 terima kasih  
 ありがとうございます  
 감사합니다  
 xixie  
 dadi  
 madloba  
 kam sab hamnida  
 rahmal  
 dindi