


phantomTM 2.0

Clear Lingual Direct Bonding Self-Ligating Bracket

2.0 hybrid



WHY A TRANSPARENT LINGUAL BRACKET?


- No dark shades from brackets visible through teeth
- Truly esthetic – also on lower teeth
- Self-ligating
- Direct bonding with fixed placement-guide guarantees correct position on tooth
- Active in final stage with larger arch wires
- Optional ceramic tubes on premolars
- Great patient comfort from protective polymer closing door



Incisors and Cuspids

Premolars

Ceramic tubes for premolars and molars



Bracket placing instructions



1. Put the placement-guide against the incisal edge



2. Remove the guide



3. Open door for wire insertion

GESTENCO

Smiles forever

www.gestenco.com
Gothenburg, Sweden