


Sunset strip™



Single-Sided Sunset strip™

Our Versatile single-sided Sunset strip™ are perfect for contouring and finishing composite restorations and for single tooth slenderization.



Part number	SU-SS-01	SU-SS-02	SU-SS-03	SU-SS-04
Grit	Super Fine	Fine	Medium	Coarse
Thickness	0.08	0.10	0.13	0.18

Double-Sided Sunset strip™

Excellent at contouring and finishing composite restorations and are especially well suited for interproximal reductions.



Part number	SU-DS-01	SU-DS-02	SU-DS-03	SU-DS-04
Grit	Super Fine	Fine	Medium	Coarse
Thickness	0.11	0.15	0.21	0.30



Serrated Sunset Strip™ - The serrated Sunset Strips™ are indispensable for separating contacts that have become fused with composite or cement.

Part number	SU-SAW-01
Thickness	0.05

GESTENCO
Smiles forever

Gestenco International AB, P.O. Box 240 67, SE-400 22 Gothenburg, Sweden
Phone +46 (0)31 81 00 35. Fax +46 (0)31 81 46 55, www.gestenco.com



Single-Sided Sunset Strip™

Our versatile single-sided Sunset Strips™ are perfect for contouring and finishing composite restorations and for single tooth slenderness.

Part number	SU-SS-01	SU-SS-02	SU-SS-03	SU-SS-04
Grit	Super Fine	Fine	Medium	Coarse
Thickness	0.08	0.10	0.13	0.18

Double-Sided Sunset Strip™

Excellent at contouring and finishing composite restorations and are especially well suited for interproximal reductions.

Part number	SU-DS-01	SU-DS-02	SU-DS-03	SU-DS-04
Grit	Super Fine	Fine	Medium	Coarse
Thickness	0.11	0.15	0.21	0.30

Serrated Sunset Strip™

The serrated Sunset Strips™ are indispensable for separating fused contacts and cement clean up.

Part number	SJ-SAW-01
Grit	No Diamond Coating
Thickness	0.05

● Super Fine Grit: 15 Micron ● Fine Grit: 30 Micron ● Medium Grit: 46 Micron ● Coarse Grit: 90 Micron

CE

GESTENCO
Smiles forever

Gestenco International AB, P.O. Box 240 67, SE-400 22 Gothenburg, Sweden
Phone +46 (0)31 81 00 35. Fax +46 (0)31 81 46 55, www.gestenco.com

Please visit us on

