

Springtime C™

Self-ligating ceramic bracket

- Interactive clear self-ligating bracket with a metal spring-clip
- Dual options for opening the spring-clip/ hole on top or pushing from the gingival side.
- Made of transparent ceramic material.
- Contoured appearance for maximum comfort
- Ceramic mesh base
- Available in Roth .022.



Passive Phase

Small round wires slide freely, initiating the tooth movement process as the archwire gently levels the teeth and coaxes them into alignment.



Expressive Phase

Square or rectangular wires are seated into the base of the slot without contacting the clip. Rotations are corrected and space closures are completed.



Active Phase

Rectangular archwires extend beyond the slot to fully engage the clip, providing the active control necessary for finishing uprighting of the roots and adjusting the torque

GESTENCO
Smiles forever

Gestenco International AB, P.O. Box 240 67, SE-400 22 Gothenburg, Sweden
Phone +46 (0)31 81 00 35. Fax +46 (0)31 81 46 55, www.gestenco.com

Springtime C™ Formula Roth*

Type	Description	Torque	Ang.	Roth .022
Centrals	Upper right	+12	+5	SPC-11-55
	left			SPC-21-55
Laterals	Upper right	+8	+9	SPC-12-55
	left			SPC-22-55
Cuspids w/hook	Upper right	0	+11	SPC-13-55HK
	left			SPC-23-55HK
Bicuspids w/hook	Upper right	-7	0	SPC-14-55HK
	left			SPC-24-55HK
Anteriors	Lower right	-1	0	SPC-40-55
	left			SPC-30-55
Cuspids w/hook	Lower right	-11	+5	SPC-43-55HK
	left			SPC-33-55HK
1st Bicuspids w/hook	Lower right	-17	0	SPC-44-55HK
	left			SPC-34-55HK
2nd Bicuspids w/hook	Lower right	-22	0	SPC-45-55HK
	left			SPC-35-55HK

GESTENCO
Smiles forever

Gestenco International AB, P.O. Box 240 67, SE-400 22 Gothenburg, Sweden
Phone +46 (0)31 81 00 35. Fax +46 (0)31 81 46 55, www.gestenco.com