Environmental samples

small samples series and low sample volumes



- IVL Swedish Environmental Research Institute in short
- Environmental samples.
- Iron Maiden.
- Small sample volumes.
- Small sample series.
- Examples of method development.





About IVL

- Independent status, owned by a foundation
- Jointly established by the Swedish State and the business sector 1966
- Operations span the whole field of sustainability
- Both research and commission services
- Alliances and networks
- Laboratories and pilot plants
- Offices in Stockholm, Gothenburg, Malmö and China
- Subsidiaries: EPD International, BASTA and IVL China



Environmental samples

















Iron Maiden ICAP 7000 series







Low sample volumes





Low sample volumes

- Filter that have been leached in water
 - Many different analyses on the same leachate
 - LOD
 - 4-5 ml sample volume.
- Percipitation samples
 - Many different analyses on the same rain sample
 - 5-10 ml sample volume



Small sample series



- Few samples.
 - 10-20 samples in one series
 - Research projects often have strained budgets
- Experimental.
- Economy.
- Practice makes perfect.
 - Learn from experiences



Diffusive Particle Samplers



- Developed at IVL
- Worldwide
- Sulfur
- Digestion
- Blanks



Filter from airsampling



- Method change
- Ion cromatograph to ICP-OES
- Kations



Thank You for your attention!

Questions?

