

Extrusion Eyelet with chains and Caplin hook

- 24K Gold plated eyelet and chain reduces possibility of allergic reactions
- Low profile limits irritation
- Mesh base assures reliable bonding

Bond the extrusion eyelet to an impacted tooth, then engage traction force.

The bonding pads come in various sizes and shapes in order to fit the anatomy of all teeth and cases. The Caplin hook is also available in Stainless Steel without chain.



Extrusion Eyelet, ovoid base



Extrusion Eyelet, round base



Extrusion Eyelet, square base



Extrusion Caplin hook



Citizen Caplin hook, Stainless Steel

Extrusion Eyelet w chain and Caplin hook	
Ovoid base w chain	EE-00-01
Round base w chain	EE-00-02
Square base w chain	EE-00-03
Caplin hook w chain	EE-00-04
Citizen Caplin hook	CIB-00-03