

# Diabetes care





# Diabetes care





## Soft Blood Glucose lancets

Soft Blood Glucose lancets are intended for blood glucose measurement. The soft lancets are universal and compatible with most lancing devices especially our Compact Eject. The range includes standard soft lancets and siliconized soft lancets available in different sizes.

### GOOD BLOOD FLOW

Triple ground needle point.

### VERSATILE

Fit most commonly used lancing devices including the Vitrex Compact Eject lancing device.

### PATIENT SAFETY:

Single-use only. Sterilised by gamma irradiation.

### CONVENIENT

Boxes have dispensing feature.

### NEW SILICONIZED NEEDLE

For less pain and smoother penetration.

### Technical Specifications

Colour Code	Specification	Pcs./Box	Article No.
Light blue	23G (0.64 mm)	200	351213
Red	28G (0.36 mm)	200	352213
Green	30G (0.30 mm)	200	354213
Dark blue	32G (0.26 mm)	200	357213

Shelf life 5 years

### Technical Specifications Siliconized Needles

Colour Code	Specification	Pcs./Box	Article No.
Red	28G (0.36 mm)	200	352413
Green	30G (0.30 mm)	200	354413

Shelf life 5 years



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 200 pcs.  
 Secondary packaging: 2,000 pcs.  
 Transit packaging: 20,000 pcs.  
 Pallet: 480,000 pcs.



## Compact Eject

The Compact Eject is a lancing device with a push button ejection. The device is clear adjustable with clear fixed end caps designed for Alternate Site Testing (AST).

### VERSATILE

Suitable for use with our Soft Lancets and other most commonly used lancet brands.

### VARIABLE PENETRATION

Nine depth settings from 0.5 to 2.1 mm

### CONVENIENT

Simple and quick activation, reactivation and non-touch lancet removal.

### Technical Specifications

Product	Description	Packaging	Article No.
Compact Eject	9-step adjustable depth	Clear plastic sleeve only	369613

Shelf life 5 years





# Medoject® Fine Pen Needles

Medoject® Fine Pen Needles are made of polypropylene hub and stainless steel tube. They are intended for Insulin application via most Insulin pens. The needles are sterile and individually packed.

## COMFORTABLE INJECTION

European, premium needle quality, triple beveled tip and precision lancet point grinding.

## THIN WALL TECHNOLOGY

Ensures optimal insulin flow.

## SILICONIZED NEEDLE

Minimizes irritation and gives a smooth penetration.

## SHORT NEEDLES

Reduce risk of a painful injection into the muscle.

## VERSATILE

Fit most commonly used insulin pens.

## COLOUR-CODING

Simple identification of different needle lengths.

## FDA CLEARED

## Technical Specifications

Colour Code	Specification	Article No.
Yellow	29G x 10 mm	CHPN2910
Pink	29G x 12.7 mm	CHPN29127
Dark blue	30G x 6 mm	CHPN3006
Light blue	30G x 8 mm	CHPN3008
Green	31G x 4 mm	CHPN3104
Dark violet	31G x 5 mm	CHPN3105
Dark blue	31G x 6 mm	CHPN3106
Light blue	31G x 8 mm	CHPN3108
Green	32G x 4 mm	CHPN3204
Dark violet	32G x 5 mm	CHPN3205
Dark blue	32G x 6 mm	CHPN3206
Green	33G x 4 mm	CHPN3304 *

Medoject Fine Pen Needles are available in many sizes.

Shelf life 5 years

For more information Contact Sales Department.

\* Not FDA cleared.



Transit packaging

Pallet

### Packaging:

Primary packaging:

100 pcs.

Transit packaging:

2,000 pcs.

Pallet:

160,000 pcs.



# Safety Pen Needles

## With Dual Automatic Protective Shields

Safety Pen Needles are intended for a safe and simple insulin injection via most insulin pens. The application is the same as normal insulin pens, but they offer an automatic double protection without extra operations after using. The coloured indicator shields in the front and back-end provide a visual, physical and audible confirmation of safety. The needles are sterile and individually packed.

### SAFETY FUNCTION

Coloured automatic safety shields on front and back-end providing double protection after using.

### COMFORTABLE INJECTION

Triple sharpened needle and top lubrication technology, which offers minimum penetration force.

### SIMPLE AND SAFE APPLICATION

Easy to trigger the safety function after injection. Same application as normal pen needles, but with a safety function both blocking top and bottom needle after use.

### VERSATILE

Uniform hub screw which fits most commonly used insulin pens.

## Technical Specifications

Colour Code	Specification	Article No.
Purple ●	30G x 5 mm	531113
Dark blue ●	30G x 8 mm	531213
Yellow ●	31G x 5 mm	531313
Orange ●	31G x 6 mm	531413
Green ●	32G x 4 mm	531513

Shelf life 5 years



**Packaging:**  
 Primary packaging: 100 pcs.  
 Transit packaging: 800 pcs.  
 Pallet: 57,600 pcs.



# Chirana® Insulin and Tuberculin Syringes

Chirana® Insulin and Tuberculin syringes are with integrated needle for single use. The needle cap is made of latex-free rubber gasket. The syringes are individually packed in boxes of 100 pcs. They are pyrogen-free and non-toxic.

**USER CONVENIENCE:** Needle permanently fixed to syringe body.

**USER SAFETY:** Integrated needle cap.

**COMFORT:** Grinded, precision lancet point and silicone treated tip for smooth penetration.

### PATIENT SAFETY:

Single-use only. Sterilised by ETO. Dead space is minimized to reduce waste and transmission of infectious agents. The dead space does not exceed the limits (EN ISO 8537:2016, art. 5.11.1).

**FDA CLEARED \***



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Secondary packaging: 2,000 pcs.  
 Pallet: 48,000 pcs.

\* FDA clearance for Insulin syringes.

## Technical Specifications

Type	Vol. ml.	Needle	Size Metric	Article No:
Tuberculin	1	29G x 1/2"	12.7 mm	CHTUB0129
Tuberculin	1	Luer		CHTUB01
Insulin U 40	1	29G x 1/2"	12.7 mm	CHINS4129
Insulin U 100	1	27G x 1/2"	12.7 mm	CHINS0127
Insulin U 100	1	29G x 5/16"	8 mm	CHINS01298
Insulin U 100	1	29G x 1/2"	12.7 mm	CHINS0129
Insulin U 100	1	30G x 1/2"	12.7 mm	CHINS0130
Insulin U 100	1	30G x 5/16"	8 mm	CHINS01308
Insulin U 100	1	31G x 4/16"	6 mm	CHINS01316
Insulin U 100	1	31G x 5/16"	8 mm	CHINS01318
Insulin U 100	1	Luer		CHINS01
Insulin U 100	0.5	29G x 1/2"	12.7 mm	CHINS00529
Insulin U 100	0.5	30G x 1/2"	12.7 mm	CHINS00530
Insulin U 100	0.5	30G x 5/16"	8 mm	CHINS005308
Insulin U 100	0.5	31G x 5/16"	8 mm	CHINS005318
Insulin U 100	0.5	31G x 4/16"	6 mm	CHINS005316
Insulin U 100	0.3	29G x 1/2"	12.7 mm	CHINS00329
Insulin U 100	0.3	30G x 1/2"	12.7 mm	CHINS00330
Insulin U 100	0.3	31G x 5/16"	8 mm	CHINS003318
Insulin U 100	0.3	31G x 4/16"	6 mm	CHINS003316

Shelf life 5 years

Chirana Syringes are available in many sizes.  
 For more information Contact Sales Department.



## Alcohol Swabs

Alcohol Swabs are an easy and convenient way to cleanse the skin and help fight bacteria at the site of application. Use as an antiseptic skin cleanser before injections and blood sampling.

The size of the Disinfection Swab is 30x60 mm. Each Alcohol Swab is individually sealed and disposable for convenience and safety.

### Technical Specifications

Description	Size	Content	Identification No.	Pcs./Box	Article No.
Alcohol Swabs	30 x 60 mm	82% Ethanol	1U8K-YGJD-J103-YQRW	100	521413
Alcohol Swabs	30 x 60 mm	82% Ethanol + 0.5% Chlorhexidine	NQ1P-3GSC-U100-G3GE	100	522413

Shelf life 5 years

Alcohol Swabs are currently approved for sales in: Sweden, Finland, UK, Germany, France and Austria.

(Updated: 2025-10. Please consult our website: [vitrexmg.com](http://vitrexmg.com) for new updates in the list of approved countries).

Regulated according to: REGULATION (EU) No 528/2012 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 22 May 2012 concerning the making available on the market and use of biocidal products.



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Transit packaging: 10,000 pcs.  
 Pallet: 400,000 pcs.



**CHIRANA T. Injecta, a.s.**  
Nám. Dr. Alberta Schweitzera 194  
916 01 Stará Turá, Slovakia  
Phone: +421 32 7709 901



**Vitrex Medical A/S**  
Vasekaer 6-8  
2730 Herlev, Denmark  
Phone: +45 4494 7011



**CHIRANA T. Injecta, s.r.o.**  
Komořanská 2148  
143 00 Praha 4  
Czech Republic  
Phone: +420 241 097 904

