





About GA Health

GA Health is a global company with its headquarters in Hong Kong and manufacturing operations in China & Ireland. It has a distribution network which now stretches globally with a name built on quality, trust and integrity.

Our mission is to improve the quality of life through the design, development and manufacture of high quality, affordable medical product, for the benefit of patients and health professionals and with respect to our employees.



About Andorate®

Andorate® is an emerging brand of disposable endoscopic accessories and devices manufactured by GA Health. Our products are designed, developed and manufactured to reflect our ongoing commitment to patient safety, infection prevention and control.

Ken McCabe,
Co-founder & President



Cindy Ye,
Co-founder & CEO



Vision & Mission

Vision

Healthcare without infection.


Mission

To improve the quality of life through the design, development and manufacture of high quality, affordable medical product, for the benefit of patients and health professionals and with respect to our employees.

Values

- Integrity
- Fairness
- Quality
- Collaboration and partnership
- Customer experience
- Simplicity

Key Milestones



For Olympus [®] GI Endoscopes	P.6
For Pentax [®] GI Endoscopes	P.16
For Fujifilm [®] GI Endoscopes	P.22
Universal	
Tubing	P.31
Polyp Trap	P.32
Bite Block	P.34
Ordering Information	
For Olympus [®] GI Endoscopes	P.38
For Pentax [®] GI Endoscopes	P.40
For Fujifilm [®] GI Endoscopes	P.41
Universal	P.42

For **Olympus**[®] GI Endoscopes



For **Pentax**[®] GI Endoscopes

For **Fujifilm**[®] GI Endoscopes




Universal

andorate[®]

Brand New Every Time

For
Olympus®
GI Endoscopes

Lubricant free Valve Kit



Suction Valve
GAR004A

- To aspirate patient fluid lightly and effectively



Air/ Water Valve
GAR004B

- To deliver air during the surgery
- To clean the lens for better vision



Biopsy Valve
GAR028

- To provide the passage for biopsy instruments



Biopsy Irrigator
GAR015

- To provide the passage for biopsy instruments
- Equipped with one-way valve to provide safe irrigation when needed

Disposable Tubing



AquaPulse®
Auxiliary Water
Connector

GAR027

- Equipped with one-way valve to provide safe irrigation



AquaPulse®
Irrigation Tubing

GAR073

For Olympus® OFP and Endo Stratus™ EGA-500 Irrigation Pump

- Equipped with one-way valve to provide safe 24-hour irrigation




AquaPulse®
Irrigation Tubing

GAR083

For Olympus® OFP, Endo Stratus™ EGA-500 and Endogator™ EGP-100 Irrigation Pump

- Equipped with one-way valve to provide safe 24-hour irrigation



Vacutore®
Air/ Water Bottle Tubing

GAR025

- To provide air supply and water supply during the endoscopy procedure
- For 24-hour use



Vacutore®
CO₂ Source Tubing
with Luer Connector

GAR047

- Connection for CO₂ insufflator
- Equipped with micro filter
- For 24-hour use

Lubricant free Valve Kit

Brand New Every Time



2-piece Valve Kit
(Suction & Air/ Water Valves)

GAR004C

For Olympus® 140/160/180/190/240/
260/290 Series GI Endoscopes, STERILE

Qty per Box 50 units



3-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves)

GAR004

For Olympus® 140/160/180/190/240/
260/290 Series GI Endoscopes, STERILE

Qty per Box 50 units



4-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves and
Auxiliary Water Connector)

GAR046

For Olympus® 160/180/190/240/
260/290 Series GI Endoscopes, STERILE

Qty per Box 50 units



Biopsy Valve

GAR028

For Olympus® & Fujifilm® GI Endoscopes, NON-STERILE

Qty per Box 100 units

GAR028-S

For Olympus® & Fujifilm® GI Endoscopes, STERILE

Qty per Box 100 units



Biopsy Irrigator

GAR015

For Olympus® & Fujifilm® GI Endoscopes, STERILE

Qty per Box 100 units

Air/ Water Bottle Tubing

Brand New Every Time



Air/ Water Bottle Tubing

GAR025

- For Olympus[®] 140/160/180/190/240/260/290 Series GI Endoscopes, STERILE

Qty per Box | 10 units



Air/ Water Bottle Tubing with CO₂ Extend

GAR025CO2EXT

- For Olympus[®] 140/160/180/190/240/260/290 Series GI Endoscopes, STERILE

Qty per Box | 10 units

Irrigation System

Brand New Every Time



Auxiliary Water Connector

GAR027

- For Olympus[®] 160/180/190/240/260/290 Series GI Endoscopes, STERILE

Qty per Box | 100 units

CO₂ Sourcing System

Brand New Every Time



CO₂ Source Tubing Adapter

GAR093

- For Olympus[®] UCR CO₂ Insufflator, NON-STERILE

Qty per Box | 100 units

Air/ Water Channel Cleaning Adapter and Valve Kit

- To ensure infection prevention & control
- To greatly simplify the pre-cleaning steps



Remarks: Excludes EUS echoendoscopes that utilize a balloon channel



Air/ Water Channel Cleaning Adapter

GAR038

For Olympus® GI Endoscopes, NON-STERILE

Qty per Box 100 units



A/W Channel Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves)

GAR3804

For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, STERILE

Qty per Box 10 units



A/W Channel Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector)

GAR3846

For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, STERILE

Qty per Box 10 units

For Olympus® GI Endoscopes

andorate®

Brand New Every Time

Cleaning Adapter and Valve Kit



Suction Channel Cleaning Adapter

GAR075

For Olympus® GI Endoscopes, NON-STERILE

Qty per Box 50 units



2-piece Cleaning Adapter Kit
(Suction & Air/Water Channel
Cleaning Adapter)

GAR3875

For Olympus® GI Endoscopes, STERILE

Qty per Box 25 units



2-piece Cleaning Adapter Kit
(Suction & Air/Water Channel
Cleaning Adapter)

GAR3875-NS

For Olympus® GI Endoscopes, NON-STERILE

Qty per Box 25 units



2-piece Cleaning Adapters + 3-piece Valve Kit
(Suction & Air/Water Channel Cleaning Adapter)
+ (Suction, Air/Water & Biopsy Valves)

GAR387504

For Olympus® 140/160/180/190/240/
260/290 Series GI Endoscopes, STERILE

Qty per Box 10 units




2-piece Cleaning Adapters + 4-piece Valve Kit
(Suction & Air/Water Channel Cleaning Adapter)
+ (Suction, Air/Water & Biopsy Valves and
Auxiliary Water Connector)

GAR387546

For Olympus® 160/180/190/240/
260/290 Series GI Endoscopes, STERILE

Qty per Box 10 units

For Pre-Cleaning Process



Direction for Use

- 1 Turn off the airflow regulator on the light source after initial rinsing of scope (per facility guidelines).
- 2 Remove the air/water valve from the endoscope.
- 3 Insert the ANDORATE® Air/ Water Channel Cleaning Adapter into the air/water cylinder on the control body of the endoscope.
- 4 Place the distal end of the scope into a basin of water.
- 5 On the light source, turn on the airflow regulator and set it to a maximum output of "high" .
- 6 Depress the ANDORATE® Air/ Water Channel Cleaning Adapter for 30 seconds.
- 7 Release the ANDORATE® Air/ Water Channel Cleaning Adapter.
- 8 Keep the distal end of the scope in the basin of water. Flush air for 10 seconds.
- 9 Turn off the airflow regulator on the light source.
- 10 Remove distal end of endoscope from basin of water.
- 11 Remove the ANDORATE® Air/ Water Channel Cleaning Adapter from endoscope after the precleaning step is completed.
- 12 Discard the ANDORATE® Air/ Water Channel Cleaning Adapter in accordance with accepted medical practice and applicable local, state and federal laws and regulations.



For Olympus® GI Endoscopes

andorate®

Brand New Every Time

GAR027



GAR046




GAR025



GAR387546

For
Pentax®
GI Endoscopes

Lubricant free Valve Kit



Suction Valve

GAR037A

- To aspirate patient fluid lightly and effectively



Air/ Water Valve

GAR037B

- To deliver air during the surgery
- To clean the lens for better vision



Biopsy Valve

GAR018

- To provide the passage for biopsy instruments

Disposable Tubing



Vacutore®
Air/ Water Bottle Tubing
GAR072

- To provide air supply and water supply during endoscopy procedure
- For 24-hour use



Vacutore®
Air/ Water Bottle Tubing
With Side Port
GAR076

- To provide air supply and water supply during endoscopy procedure
- For 24-hour use



Vacutore®
CO₂ Source Tubing
with Luer Connector
GAR047

- Connection for CO₂ insufflator
- Equipped with micro filter
- For 24-hour use



AquaPulse®
Auxiliary Water Connector
GAR048

- Equipped with one-way valve to provide safe irrigation



AquaPulse®
Irrigation Tubing
GAR073

For Olympus® OFP and Endo Stratus™ EGA-500 Irrigation Pump


- Equipped with one-way valve to provide safe 24-hour irrigation



AquaPulse®
Irrigation Tubing
GAR083

For Olympus® OFP, Endo Stratus™ EGA-500 and Endogator™ EGP-100 Irrigation Pump

- Equipped with one-way valve to provide safe 24-hour irrigation



Lubricant free Valve Kit

Brand New Every Time



2-piece Valve Kit
(Suction & Air/ Water Valves)

GAR037C

For Pentax® 90/ i10 Series GI Endoscopes,
STERILE

Qty per Box 50 units



3-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves)

GAR037

For Pentax® 90/ i10 Series GI Endoscopes,
STERILE

Qty per Box 50 units



4-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves and
Auxiliary Water Connector)

GAR081

For Pentax® 90/ i10 Series GI Endoscopes,
STERILE

Qty per Box 50 units



Biopsy Valve

GAR018

For Pentax® GI Endoscopes, NON-STERILE

Qty per Box 100 units

GAR018-S

For Pentax® GI Endoscopes, STERILE

Qty per Box 100 units

For **Pentax®** GI Endoscopes

andorate®

Brand New Every Time

Air/ Water Bottle Tubing



Air/ Water Bottle Tubing

GAR072

- For Pentax[®] 90 Series GI Endoscopes, STERILE

Qty per Box | 10 units



Air/ Water Bottle Tubing with CO₂ Extend

GAR072CO2EXT

- For Pentax[®] 90 Series GI Endoscopes, STERILE

Qty per Box | 10 units



Air/ Water Bottle Tubing with Side Port

GAR076

- For Pentax[®] 90 Series GI Endoscopes, STERILE

Qty per Box | 10 units

Irrigation System



Auxiliary Water Connector

GAR048

- For Pentax[®] 90/ i10 Series GI Endoscopes, STERILE

Qty per Box | 100 units



GAR081



GAR072



GAR076


For Pentax® GI Endoscopes

andorate®

Brand New Every Time

For
Fujifilm[®]
GI Endoscopes

Lubricant free Valve Kit



Suction Valve

GAR087A

- To aspirate patient fluid lightly and effectively



Air/ Water Valve

GAR087B

- To deliver air during the surgery
- To clean the lens for better vision



Biopsy Valve

GAR028

- To provide the passage for biopsy instruments



Biopsy Irrigator

GAR015

- To provide the passage for biopsy instruments
- Equipped with one-way valve to provide safe irrigation when needed

Disposable Tubing



**AquaPulse®
Auxiliary Water
Connector**

GAR064

- Equipped with one-way valve to provide safe irrigation



**AquaPulse®
Irrigation Tubing**

GAR073

For Olympus® OFP and Endo Stratus™ EGA-500 Irrigation Pump

- Equipped with one-way valve to provide safe 24-hour irrigation




**AquaPulse®
Irrigation Tubing**

GAR083

For Olympus® OFP, Endo Stratus™ EGA-500 and Endogator™ EGP-100 Irrigation Pump

- Equipped with one-way valve to provide safe 24-hour irrigation



**Vacutore®
Air/ Water Bottle Tubing**

GAR090

- To provide air supply and water supply during the endoscopy procedure
- For 24-hour use



**Vacutore®
CO2 Source Tubing
with Luer Connector**

GAR047

- Connection for CO2 insufflator
- Equipped with micro filter
- For 24-hour use

Lubricant free Valve Kit

Brand New Every Time



2-piece Valve Kit
(Suction & Air/ Water Valves)

GAR087C

For Fujifilm®700 Series GI Endoscopes,
STERILE

Qty per Box 50 units



3-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves)

GAR087

For Fujifilm®700 Series GI Endoscopes,
STERILE

Qty per Box 50 units



4-piece Valve Kit
(Suction, Air/ Water & Biopsy Valves and
Auxiliary Water Connector)

GAR092

For Fujifilm®700 Series GI Endoscopes,
STERILE

Qty per Box 50 units



Biopsy Valve

GAR028

For Fujifilm® & Olympus® GI Endoscopes, NON-STERILE

Qty per Box 100 units

GAR028-S

For Fujifilm® & Olympus® GI Endoscopes, STERILE

Qty per Box 100 units



Biopsy Irrigator

GAR015

For Fujifilm® & Olympus® GI Endoscopes, STERILE

Qty per Box 100 units

Air/ Water Bottle Tubing

Brand New Every Time



Air/ Water Bottle Tubing

GAR090

- For Fujifilm[®] 700 Series GI Endoscopes, STERILE

Qty per Box 10 units



Air/ Water Bottle Tubing with CO₂ Extend

GAR090CO2EXT

- For Fujifilm[®] 700 Series GI Endoscopes, STERILE

Qty per Box 10 units



Air/ Water Bottle Tubing

GAR088

- For Fujifilm[®] 500/ 600 Series GI Endoscopes, STERILE

Qty per Box 10 units



Air/ Water Bottle Tubing with CO₂ Extend

GAR088CO2EXT

- For Fujifilm[®] 500/ 600 Series GI Endoscopes, STERILE

Qty per Box 10 units

Irrigation System

Brand New Every Time



Auxiliary Water Connector

GAR064

- For Fujifilm[®] 500/ 600/ 700 Series GI Endoscopes, STERILE

Qty per Box 100 units

CO₂ Sourcing System

Brand New Every Time



NEW

CO₂ Source Tubing Adapter

GAR0108

- For Fujifilm[®] GW-100 Insufflator, NON-STERILE

Qty per Box 100 units

Air/ Water Channel Cleaning Adapter and Valve Kit

- To ensure infection prevention & control
- To greatly simplify the pre-cleaning steps

NEW



Remarks: Excludes EUS echoendoscopes that utilize a balloon channel



Air/ Water Channel Cleaning Adapter

GAR105

For Fujifilm® 700 Series GI Endoscopes,
NON-STERILE

Qty per Box 100 units



A/W Channel Cleaning Adapter
+ 3-piece Valve Kit
(Suction, Air/Water & Biopsy Valves)

GAR105087

For Fujifilm® 700 Series GI Endoscopes,
STERILE

Qty per Box 10 units



A/W Channel Cleaning Adapter
+ 4-piece Valve Kit
(Suction, Air/Water & Biopsy Valves
and Auxiliary Water Connector)

GAR105092

For Fujifilm® 700 Series GI Endoscopes,
STERILE


Qty per Box 10 units

For Fujifilm® GI Endoscopes

andorate®

Brand New Every Time

For Pre-Cleaning Process



Direction for Use

- 1 Turn off the airflow regulator on the light source after initial rinsing of scope (per facility guidelines).
- 2 Remove the air/water valve from the endoscope.
- 3 Insert the ANDORATE® Air/ Water Channel Cleaning Adapter into the air/water cylinder on the control body of the endoscope.
- 4 Place the distal end of the scope into a basin of water.
- 5 On the light source, turn on the airflow regulator and set it to a maximum output of "high" .
- 6 Depress the ANDORATE® Air/ Water Channel Cleaning Adapter for 30 seconds.
- 7 Release the ANDORATE® Air/ Water Channel Cleaning Adapter.
- 8 Keep the distal end of the scope in the basin of water. Flush air for 10 seconds.
- 9 Turn off the airflow regulator on the light source.
- 10 Remove distal end of endoscope from basin of water.
- 11 Remove the ANDORATE® Air/ Water Channel Cleaning Adapter from endoscope after the precleaning step is completed.
- 12 Discard the ANDORATE® Air/ Water Channel Cleaning Adapter in accordance with accepted medical practice and applicable local, state and federal laws and regulations.



For Fujifilm® GI Endoscopes

andorate®

Brand New Every Time

The image features a minimalist abstract design. A central, bright orange, rounded, leaf-like shape is positioned horizontally. It is flanked by two large, light grey, rounded shapes that appear to be part of a larger, stylized letter 'U' or a similar abstract form. The word "Universal" is written in a clean, white, sans-serif font, centered within the orange shape.

Universal

Air/Water Bottle Tubing

Brand New Every Time



CO₂ Source Tubing with Luer Connector

GAR047
NON-STERILE

Qty per Box 10 units



CO₂ Extend Tubing with Luer connector

GAR091
NON-STERILE

Qty per Box 20 units

Irrigation Tubing

Brand New Every Time



Irrigation Tubing

GAR073
For Olympus[®] OFP and
Endo Stratus[™] EGA-500 Irrigation Pump, STERILE

Qty per Box 10 units




Irrigation Tubing

GAR083
For Olympus[®] OFP, Endo Stratus[™] EGA-500 and
Endogator[™] EGP-100 Irrigation Pump, STERILE

Qty per Box 10 units

Click4Catch™ Click into Place Polyp Trap with Four Removable Capture Chambers

- Ergonomic Design
- Accurate specimen identification
- Minimize the exposure to collected biomaterials



andorate® Polyp Trap

Brand New Every Time



Click4Catch™ Click into Place Polyp Trap with Four Removable Capture Chambers

GAR082

NON-STERILE

Qty per Box 12 units



ThomasTrap™ Polyp Trap with Two Removable Chambers

GAR040

NON-STERILE

Qty per Box 20 units



Four Chamber Polyp Trap

GAR023

NON-STERILE

Qty per Box 12 units



PearlCatch™ Single Chamber Polyp Trap

GAR014

NON-STERILE

Qty per Box 50 units


Universal - Polyp Trap

andorate®

Brand New Every Time

Bite Block

- One size fits all endoscopes



andorate® Bite Block

Brand New Every Time



Bite Block 66Fr with Strap

SDPC0210
NON-STERILE

Qty per Box 100 units



Bite Block 66Fr with Strap and Tooth Protector

SDPC0210F
NON-STERILE


Qty per Box 100 units



Bite Block 55.5Fr with Strap

SDPC0209A
NON-STERILE

Qty per Box 100 units



Bite Block 55.5Fr with Strap and Tooth Protector

SDPC0209
NON-STERILE


Qty per Box 100 units



Bite Block 66Fr with Oxygen Connection

GAR044
NON-STERILE

Qty per Box 100 units



Bite Block 66Fr with Oxygen Connection and Tooth Protector

GAR044F
NON-STERILE

Qty per Box 100 units





Ordering Information






Ordering Information

For Olympus® GI Endoscopes

Product Code	Description	Qty Per Box
GAR004	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR004A	Suction Valve For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	100 units
GAR004B	Air/ Water Valve For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	100 units
GAR004C	2-piece Valve Kit (Suction & Air/Water Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR015	Biopsy Irrigator For Olympus® and Fujifilm® GI Endoscopes, Sterile	100 units
GAR025	Vacutore® Air/ Water Bottle Tubing For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR025CO2EXT	Vacutore® Air/ Water Bottle Tubing with CO2 Extend For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR027	Aquapulse® Auxiliary Water Connector For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	100 units
GAR028	Biopsy Valve For Olympus® and Fujifilm® GI Endoscopes, Non-Sterile	100 units
GAR028-S	Biopsy Valve For Olympus® and Fujifilm® GI Endoscopes, Sterile	100 units
GAR038	Air/ Water Channel Cleaning Adapter For Olympus® GI Endoscopes, Non-Sterile	100 units
GAR3804	A/W Channel Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR3846	A/W Channel Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR3875	2-piece Cleaning Adapter Kit (Suction & Air/Water Channel Cleaning Adapter) For Olympus® GI Endoscopes, Sterile	25 units
GAR3875-NS	2-piece Cleaning Adapter Kit (Suction & Air/Water Channel Cleaning Adapter) For Olympus® GI Endoscopes, Non-Sterile	25 units
GAR387504	2-piece Cleaning Adapters + 3-piece Valve Kit (Suction & Air/Water Channel Cleaning Adapter) + (Suction, Air/Water & Biopsy Valves) For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR387546	2-piece Cleaning Adapters + 4-piece Valve Kit (Suction & Air/Water Channel Cleaning Adapter) + (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	10 units
GAR046	4-piece Valve Kit (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector) For Olympus® 160/180/190/240/260/290 Series GI Endoscopes, Sterile	50 units
GAR075	Suction Channel Cleaning Adapter For Olympus® GI Endoscopes, Non-Sterile	50 units
GAR093	Vacutore® CO2 Source Tubing Adapter For Olympus® UCR CO2 Insufflator, Non-Sterile	100 units

Ordering Information

For **Olympus®** GI Endoscopes

COMBINATION	GAR004C 2-piece Valve Kit (Suction & Air/ Water Valves)	GAR027 Aquapulse® Auxiliary Water Connector	GAR028-S Biopsy Valve	GAR038 Air/ Water Channel Cleaning Adapter	GAR075 Suction Channel Cleaning Adapter
					
GAR004	●		●		
GAR046	●	●	●		
GAR3804	●		●	●	
GAR3846	●	●	●	●	
GAR3875				●	●
GAR387504	●		●	●	●
GAR387546	●	●	●	●	●

For Pentax® GI Endoscopes

Product Code	Description	Qty Per Box
GAR018	Biopsy Valve For Pentax® GI Endoscopes, Non-Sterile	100 units
GAR018-S	Biopsy Valve For Pentax® GI Endoscopes, Sterile	100 units
GAR037	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Pentax® 90/ i10 Series GI Endoscopes, Sterile	50 units
GAR037A	Suction Valve For Pentax® 90/ i10 Series GI Endoscopes, Sterile	100 units
GAR037B	Air/ Water Valve For Pentax® 90/ i10 Series GI Endoscopes, Sterile	100 units
GAR037C	2-piece Valve Kit (Suction & Air/Water Valves) For Pentax® 90/ i10 Series GI Endoscopes, Sterile	50 units
GAR048	Aquapulse® Auxiliary Water Connector For Pentax® 90/ i10 Series GI Endoscopes, Sterile	100 units
GAR072	Vacutore® Air/ Water Bottle Tubing For Pentax® 90 Series GI Endoscopes, Sterile	10 units
GAR072CO2EXT	Vacutore® Air/ Water Bottle Tubing with CO2 Extend For Pentax® 90 Series GI Endoscopes, Sterile	10 units
GAR076	Vacutore® Air/ Water Bottle Tubing with Side Port For Pentax® 90 Series GI Endoscopes, Sterile	10 units
GAR081	4-piece Valve Kit (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector) For Pentax® 90/ i10 Series GI Endoscopes, Sterile	50 units

For Fujifilm® GI Endoscopes

Product Code	Description	Qty Per Box
GAR015	Biopsy Irrigator For Fujifilm® & Olympus® GI Endoscopes, Sterile	100 units
GAR028	Biopsy Valve For Fujifilm® & Olympus® GI Endoscopes, Non-Sterile	100 units
GAR028-S	Biopsy Valve For Fujifilm® & Olympus® GI Endoscopes, Sterile	100 units
GAR064	Aquapulse® Auxiliary Water Connector For Fujifilm® 500/ 600/ 700 Series GI Endoscopes, Sterile	100 units
GAR087	3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Fujifilm®700 Series GI Endoscopes, Sterile	50 units
GAR087A	Suction Valve For Fujifilm®700 Series GI Endoscopes, Sterile	100 units
GAR087B	Air/ Water Valve For Fujifilm®700 Series GI Endoscopes, Sterile	100 units
GAR087C	2-piece Valve Kit (Suction & Air/Water Valves) For Fujifilm®700 Series GI Endoscopes, Sterile	50 units
GAR088	Vacutore® Air/ Water Bottle Tubing For Fujifilm® 500/ 600 Series GI Endoscopes, Sterile	10 units
GAR088CO2EXT	Vacutore® Air/ Water Bottle Tubing with CO2 Extend For Fujifilm® 500/ 600 Series GI Endoscopes, Sterile	10 units
GAR090	Vacutore® Air/ Water Bottle Tubing For Fujifilm® 700 Series GI Endoscopes, Sterile	10 units
GAR090CO2EXT	Vacutore® Air/ Water Bottle Tubing with CO2 Extend For Fujifilm® 700 Series GI Endoscopes, Sterile	10 units
GAR092	4-piece Valve Kit (Suction, Air/Water & Biopsy Valves & Auxiliary Water Connector) For Fujifilm®700 Series GI Endoscopes, Sterile	50 units
● GAR105	Air/ Water Channel Cleaning Adapter For Fujifilm®700 Series GI Endoscopes, Non-Sterile	100 units
● GAR105087	A/W Channel Cleaning Adapter + 3-piece Valve Kit (Suction, Air/Water & Biopsy Valves) For Fujifilm®700 Series GI Endoscopes, Sterile	10 units
● GAR105092	A/W Channel Cleaning Adapter + 4-piece Valve Kit (Suction, Air/Water & Biopsy Valves and Auxiliary Water Connector) For Fujifilm®700 Series GI Endoscopes, Sterile	10 units
● GAR108	Vacutore® CO2 Source Tubing Adapter For Fujifilm® GW-100 Insufflator, Non-Sterile	100 units

- New/ Advanced SKUs

Ordering Information

For **Fujifilm®** GI Endoscopes

COMBINATION	GAR087C 2-piece Valve Kit (Suction & Air/ Water Valves)	GAR028-S Biopsy Valve	GAR064 Aquapulse® Auxiliary Water Connector	GAR105 Air/ Water Channel Cleaning Adapter
				
GAR087	●	●		
GAR092	●	●	●	
GAR105087	●	●		●
GAR105092	●	●	●	●

Universal - Tubing

Product Code	Description	Qty Per Box
GAR047	Vacutore® CO2 Source Tubing with Luer Connector Non-Sterile	10 units
GAR073	AquaPulse® Irrigation Tubing For Olympus® OFP and Endo Stratus™ EGA-500 Irrigation Pump, Sterile	10 units
GAR083	AquaPulse® Irrigation Tubing For Olympus® OFP, Endo Stratus™ EGA-500 and Endogator™ EGP-100 Irrigation Pump, Sterile	10 units
GAR091	Vacutore® CO2 Entend Tubing Luer Connector Non-Sterile	20 units

andorate®

Brand New Every Time

Universal - Polyp Trap

Product Code	Description	Qty Per Box
GAR014	PearlCatch™ Single Chamber Polyp Trap Non-Sterile	50 units
GAR023	Four Chamber Polyp Trap Non-Sterile	12 units
GAR040	ThomasTrap™ Polyp Trap with Two Removable Chambers Non-Sterile	20 units
GAR082	Click4Catch™ Click into Place Polyp Trap with Four Removable Capture Chambers Non-Sterile	12 units

Universal - Bite Block

GAR044	Bite Block 66Fr with Oxygen Connection Non-Sterile	100 units
GAR044F	Bite Block 66Fr with Oxygen Connection and Tooth Protector Non-Sterile	100 units
SDPC0209	Bite Block 55.5Fr with Strap and Tooth Protector Non-Sterile	100 units
SDPC0209A	Bite Block 55.5Fr with Strap Non-Sterile	100 units
SDPC0210	Bite Block 66Fr with Strap Non-Sterile	100 units
SDPC0210F	Bite Block 66Fr with Strap and Tooth Protector Non-Sterile	100 units



Manufacturer


GA Health Company Limited

Unit 18, 21/F, Metropole Square, 2 On Yiu Street, Shatin, N.T., Hong Kong
Tel: +852 2833 9010 info@GAhealth.net www.GAhealth.net

European Headquarters

GA Health Medical Devices Limited

Forest Park Business Campus, Mullingar, N91 VEOH, Co. Westmeath, Ireland
info@GAhealth.net www.GAhealth.net



All trademarks are the property of their respective owners.

GAH-SO-MKT-C103 (ANDORATE Ver1.0)