

TROPApuck
1 g Molecular sieve
(Formulation M-3003-398)



Applications

Fields of use

- Pharmaceutical, Nutraceutical, Diagnostics, Medical Devices, Automotive and Electronics
- Industrial Applications

Options

- Specific packaging upon request
- Standard colour is natural white – other colours upon request
- Formulation can be customized

Benefits

Puck Design

- Can be easily integrated in an existing machine
- No loose desiccant
- Easy processing – Ideal also for high-speed automations

Activ-Polymer Technology

- 100% utilization of the desiccant material capabilities
- High absorption capacity
- Maintains a low relative humidity
- Residual moisture in dosage can be removed
- Extends product's shelf-life

Regulatory and Quality

- ISO-9001 certified company
- Listed in the DMF
- GMP manufacturing upon request
- Complies with Regulation (EU) No 10/2011

Characteristics

Please note that these values are average measurements

Puck

Diameter	20,0 mm
Thickness	4,0 mm
Weight	1,5 grams
Colour	White

Activ-Polymer Technology

Desiccant Type	4A Molecular Sieve
Minimum Capacity	200 mg

ART. # 03 02 0038 MS20 WH

In Colour Clear – 4.960 Pieces per Box – 12 Boxes – 59.520 Pieces per Pallet