

Selenozyme AB

559250-2271


Nanna Svartz väg 2

A Working Lab - Campus Solna

www.selenozyme.com

Product Specification

Product No.	HTRXR2-100 (-200, -500, -1000)
Lot No.	HTRXR2200701
Product Description:	Recombinant Human Thioredoxin Reductase 2 (Non-tagged; No signal peptide; 488 residues) 1 mg/mL, 100 (200, 500, 1000) μ L, 15 U per mg enzyme, in TE buffer with 50% glycerol
Specifications:	1 U of enzyme reduces 1 μ mol DTNB per min in 0.5 mL standard DTNB assay containing 2.5 mM DTNB and 0.3 mM NADPH in TE buffer (50 mM Tris-HCl, 2 mM EDTA, pH 7.5)
References:	Cheng Q, Arnér ES. (2017) Selenocysteine Insertion at a Predefined UAG Codon in a Release Factor 1 (RF1)-depleted Escherichia coli Host Strain Bypasses Species Barriers in Recombinant Selenoprotein Translation. J Biol Chem. 292:5476-5487. Arnér ES, Holmgren A. (2001) Measurement of thioredoxin and thioredoxin reductase. Curr Protoc Toxicol. Chapter 7:Unit 7.4.
Storage:	Shipped in cold pack by courier service. Store at -20° C upon delivery.
Examples of Use:	DTNB assay; Insulin-coupled Thioredoxin assay; Inhibitor screening; etc.
NOTE!	Centrifuge the tube briefly before opening the lid
Purity and Gel Analysis:	Purity is \geq 95% as judged by gel analysis



Coomassie stained SDS PAGE analysis of 10 μ g HTRXR2-100 product

Lane 1: Molecular weight marker (SeeBlue™ Plus2)

Lane 2: HTRXR2-100

